CiD
Carbostent® Implantable Devices

PERIPHERAL VASCULAR

CFRE  **BTK™**

AMPHILIMUS™ ELUTING PERIPHERAL STENT

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of  **Alvimedica** group
www.alvimedica.com

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CRE8™ BTK

Infrapopliteal Sirolimus-eluting stent coated with α Carbofilm™ and rapid-exchange balloon catheter

1. DESCRIPTION

The CRE8™ BTK device consists of a Sirolimus-eluting infrapopliteal stent firmly held on the distal end of a semicompliant balloon catheter.

The stent is a flexible implantable device that can be expanded using a PTA catheter. The stent is made of cobalt-chromium alloy (L605) and is coated with α Carbofilm™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs. Coating the substrate with α Carbofilm™ provides it with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

The outer surface of the stent has dedicated grooves, fully coated with α Carbofilm™, for containing the pharmaceutical formulation Amphilimus™, which is composed of the drug Sirolimus and a mixture of long-chain fatty acids.

The specific drug dosage for the vessel surface treated is approximately 0.9 μ g/mm², which corresponds to a minimum dose of 50 μ g on the smaller stent (2.25x8mm) and a maximum dose of 395 μ g on the larger stent (4.0x38mm - 3.5x46mm).

Two radio-opaque platinum markers at either end of the stent enable accurate positioning on the lesion to be treated.

The rapid-exchange type balloon catheter provides a safe means for delivering the infrapopliteal stent to the lesion to be treated.

The distal portion of the catheter consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radio-opaque markers, located outside the usable length of the balloon, enable accurate placement across the stenosis.

The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.

Two depth indicators placed at 90 and 100 cm from the distal end, show when the balloon comes out of the guiding catheter.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

The Manufacturer directly manufactures the CRE8™ BTK device and carries out all quality controls both during production and on the finished product, according to the norms of Good Manufacturing Practice.

2. INTENDED USE

The stent is indicated to help restore patency to peripheral vessels (e.g., infrapopliteal vessels such as the anterior tibial artery, posterior tibial artery, peroneal artery, tibio-peroneal trunk) in patients suffering from symptomatic chronic ischemia.

3. INDICATIONS

The stent is indicated in the following cases:

- Treatment of patients suffering from symptomatic chronic ischemia at infrapopliteal level, characterized by pain of ischemic origin, claudication, ulcers and/or trophic lesions with risk of lower limb amputation, in order to improve or resolve the above problems.
- Treatment of obstructive lesions of the infrapopliteal vessels in patients meeting the requirements for a PTA procedure.
- Primary peripheral vessel angioplasty with elective indication of stent implantation in vessels whose diameter is compatible with the models in the catalog.
- Short- and long-term unsatisfactory results arising from PTA procedures

4. CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations:

- Lesions considered untreatable with PTA or other interventional techniques.
- Pregnant women.
- Disorders / allergies that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to the contrast medium or drugs used during the procedure.
- Stenosis of an unprotected vessel.
- Resistant lesions that cannot be dilated.
- Lesions on a vessel with reference diameter <2.25 mm.
- Aneurysms immediately adjacent to the stent implantation site.
- Patients with known hypersensitivity or allergies to Sirolimus, fatty acids (such as stearic acid, palmitic acid, behenic acid) or the metal components of the stent.

5. MODEL

Each CRE8™ BTK device is identified by a model code and a lot number; available product codes are listed in Table 1. The code is composed of the letters IC, followed by two letters identifying the type of device, two or three digits indicating the nominal diameter of expansion of the stent and a further two digits indicating the length of the stent. The batch number allows trackability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate trackability of the device at the user end, the product code is printed on the adhesive labels on each box; these labels can be attached to the implant patient's medical chart.

6. HOW SUPPLIED

The CRE8™ BTK device is supplied sterile in a foil pouch that must not be brought into a sterile area.

N.B.: The foil pouch is the only sterile barrier.

The manufacturer performs sterilization using a mixture of ethylene oxide and CO₂. Sterility and chemical stability are guaranteed while the package remains intact and correctly stored, until the expiry date printed on the packaging (USE BEFORE DATE).

7. STORAGE

Store in a cool, dry place at a temperature of around 25 °C; excursions permitted to 30 °C.

8. WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause contamination risks with the device and cause infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The CRE8™ BTK device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce alterations.
- **Do not manipulate, touch, or handle the stent with your fingers**, which may cause damage, contamination, or stent dislodgement from the delivery balloon.
- **The stent should not come into contact with liquids prior to preparation and positioning. However, if it is absolutely necessary to irrigate the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).**
- The CRE8™ BTK device is designed to perform as a system. Do not use its components separately.
- The CRE8™ BTK device is indicated for use in combination with PTA. The delivery catheter is not a peripheral dilatation catheter; it should be used for stent deployment only.
- Do not use the CRE8™ BTK device if its packaging has been opened or damaged or not correctly stored, or if it is past its expiry date. In such cases, its sterility is not guaranteed.
- Do not use the CRE8™ BTK if, during handling, due to torsion or resistance to insertion, the proximal portion of the catheter has kinks or bends; in such cases, do not try to straighten the catheter.
- The CRE8™ BTK device should be used only by physicians trained specifically to perform percutaneous transluminal angioplasty (PTA) and peripheral stent implantation.
- A surgery team must be available for a possible intervention.
- The CRE8™ BTK should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- The physician in charge may decide to perform pre-dilatation for stent deployment in lesions that allow such an approach.
- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- Procedural success, in terms of restoration of peripheral vascularity, may depend not only on the effectiveness of the device, but also on the characteristics of the treated vessel segment and the downstream perfusion capacity of the district.
- The stent should be implanted in the target lesion using its delivery system.
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
- Do not attempt to clean or re-sterilise devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- Follow the manufacturer's instructions when using accessories (introducer, guiding catheter, guide wire, haemostatic valve).
- If unusual resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the introducer, the guiding catheter and the complete CRE8™ BTK device as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.
- If the patient has multiple lesions in a single vessel, it is advisable to first treat the distal lesion, then the proximal lesions. This order will reduce the need to cross the proximal stent during implantation of the distal stent and therefore reduces the risk of displacing the proximal stent.
- Implantation of an infrapopliteal stent may cause dissection of the vessel distal and/or proximal to the stent and may also cause sudden occlusion of the vessel, making an additional intervention necessary (surgery, further dilatation, placement of additional stents, or other procedures).
- If the stent is lost in the vessel, recovery procedures can be initiated. The procedures can, however, cause injury to vessels and/or to the vessel access site.
- During catheter inflation do not exceed the Rated Burst Pressure.
- Never use air or any other gas to inflate the balloon.
- The physician should inform the patient that if the course of the disease makes the use of bandages and/or physiotherapy necessary, he/she must avoid any local stress to the stent implantation area.

9. PRECAUTIONS

- Use of a gauged inflation device is strongly recommended.
- Where vessels are excessively tortuous and there are also diffuse atheromatous plaques, proximal atherosclerosis, the catheter may be difficult to advance. In such cases incorrect handling could lead to dissection or rupture of the vessel.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.
- When multiple stents are required, stent materials should be of similar composition.
- To minimise the risk of stent migration in the presence of a strong magnetic field, as a precautionary measure, patients who have received a stent should not undergo magnetic resonance imaging (MRI) until the vessel has fully healed (approx. 8 weeks). The MRI image quality may be compromised if the area of interest coincides with or is close to the position of the stent. Data from the literature¹, however, have shown no side effects such as RF-induced heating and displacement for systems operating at 1.5 Tesla or less, as a result of MRI on devices made of Cr-Co alloys.

10. DRUG REGIME

Studies in clinical literature indicate the need to administer anticoagulant therapy during the procedure and anti-platelet treatment after the procedure.

¹ Shellock F.G., Shellock V. J. "Metallic Stents: Evaluation of MR Imaging Safety", Am J Roentgenol 1999; 173:543-547

Antiplatelet treatment during the period after the procedure involves the administration of:

- Thienopyridines (Ticlopidine or Clopidogrel or Prasugrel in the doses indicated in the respective instruction booklets or according to medical prescription)
- Aspirin at a minimum amount of 100 mg/day indefinitely

11. INTERACTIONS WITH DRUGS

Although no specific clinical data are available, certain drugs, such as Tacrolimus, which act through the same binding protein (FKBP), may interfere with the efficacy of Sirolimus. There have been no studies on the interaction between drugs. Sirolimus is metabolized by CYP3A4. Strong CYP3A4 inhibitors may cause greater exposure to Sirolimus up to levels associated with systemic effects, especially in the case of multiple implant. Systemic exposure to Sirolimus if the patient is taking concomitant systemic immunosuppressive drugs must also be considered.

12. DIRECTIONS FOR USE

Stent + delivery catheter preparation

Before commencing the angioplasty procedure, administer heparin and check that the patient's Activated Clotting Time (ACT) is greater than 300 seconds.

During preparation of the system, take special care not to let any liquid come into contact with the stent. However, if it is absolutely necessary to irrigate the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).

The stent should be implanted in the target lesion using its delivery system.

After having examined the package for any damage, remove the CRE8™ BTK device and bring it into a sterile area.

- Check that the delivery catheter does not present kinks, bends or other damage.
- Carefully slide the protective cover off the stent by grasping the cover at the distal end. Check that the stent is intact and well centred on the balloon.
Damage to the CRE8™ BTK device may impair its performance.
- Flush the guidewire lumen with heparin/saline mixture. **Warning: do not place the stent into contact with liquids while flushing the guidewire lumen.**
- Prepare the inflation device according to the manufacturer's instructions.
- Purge any air from the balloon the stent is mounted on as follows:
 - Fill the inflation device with 4 ml of contrast medium;
 - After connecting the inflation device to the delivery catheter Luer connector, point the distal tip of the catheter downwards;
 - Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to normal** as the system fills with contrast medium.
 - Without introducing air, repeat step 3 and aspirate for 10-15 seconds until bubbles no longer appear.

Stent insertion

Introducers with a 4F diameter (internal diameter 1.35 mm) or larger are compatible for use with the CRE8™ BTK device.

Guiding catheters with a 5F diameter (internal diameter 1.47 mm) or larger are compatible for use with the CRE8™ BTK device.

Guidewires with a diameter of 0.014 inches (0.365 mm) or less are compatible for use with the CRE8™ BTK device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- Keep the delivery catheter the stent is pre-mounted on at ambient pressure during insertion.
- Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- Ensure that the haemostasis valve is completely open before introducing the system.
- Slowly advance the system over the guidewire to allow retrograde filling of the guide catheter to the blood; keep advancing the system until the stent reaches the treatment site.
- This should be done by observing, under fluoroscopy, the position of the balloon markers in relation to the stent markers.

If unusual resistance is encountered at any time, do not force the system: withdraw the introducer (or the guiding catheter) and the complete CRE8™ BTK device as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.

Stent implantation and expansion

- Position the stent and, with the help of the balloon's radiopaque markers, make sure it is correctly positioned in relation to the lesion to be treated and completely covered.
- Expand the stent by slowly inflating the balloon to the nominal diameter. Table 2 lists the diameters of the balloon as inflation pressure varies (column I), for the seven nominal diameter groups: 2.25 mm (column II), 2.5 mm (column III), 2.75 mm (column IV), 3.0 mm (column V), 3.5 mm (column VI), 4.0 mm (column VII) and 4.5 mm (column VIII). Values on dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a minimum balloon inflation pressure of about 5 atm.
- Stent elastic recoil is between 2 and 7%* depending on stent model and expansion diameter.**
* NOTE: These values are obtained from the results of *in-vitro* testing
- Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation.**
- Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTA procedure.
- If post-dilatation of the stent is required, a PTA balloon can be used. Use extreme care:
 - not to expand the 2.25 mm nominal diameter stent over 2.55 mm.
 - not to expand the 2.50 mm nominal diameter stent over 3.05 mm.
 - not to expand the 2.75 mm nominal diameter stent over 3.05 mm.
 - not to expand the 3.00 mm nominal diameter stent over 3.85 mm.
 - not to expand the 3.50 mm nominal diameter stent over 3.85 mm.
 - not to expand the 4.00 mm nominal diameter stent over 5.05 mm.
 - not to expand the 4.50 mm nominal diameter stent over 5.05 mm.

Optimal expansion requires the stent to be in full contact with the arterial wall,

**so that the stent diameter equals the diameter of the reference vessel.
ENSURE THAT THE STENT IS FULLY EXPANDED.**

13. POTENTIAL ADVERSE EVENTS

Implantation of a stent may be associated with the following adverse events:

- Death
- Dissection, perforation, rupture and/or damage to the blood vessel
- Access-site haematoma
- Distal and/or stent embolism
- Access-site bleeding
- Access-site infection and pain
- Stent occlusion
- Acute occlusion of the vessel
- Chronic occlusion of the vessel (restenosis of the lesion at stent level)
- Allergic reaction (to the contrast medium, to the drugs used during the procedure or to the materials used to make the stent).
- Thrombosis (acute, subacute or late)
- Total or partial amputation of the limb
- Arteriovenous fistula
- Femoral pseudoaneurysm
- MRI Artefacts
- Systemic bleeding

14. LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. Under such circumstances, the customer shall deliver the device to The Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty consists exclusively in the replacement of the defective device with another device of the same type or equivalent, from the same Manufacturer.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for return of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package.

Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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BLANK PAGE

BLANK PAGE

SIROLIMUS ELUTING
CARBOFILM™ COATED PERIPHERAL STENT
ON RX BALLOON CATHETER


TABLE 1		Stent length							
		8mm	12mm	16mm	20mm	25mm	31mm	38mm	46mm
Nominal stent diameter	2.25 mm	ICLK22508	ICLK22512	ICLK22516	ICLK22520	ICLK22525	ICLK22531	-	-
	2.5 mm	ICLK2508	ICLK2512	ICLK2516	ICLK2520	ICLK2525	ICLK2531	ICLK2538	ICLK2546
	2.75 mm	ICLK27508	ICLK27512	ICLK27516	ICLK27520	ICLK27525	ICLK27531	ICLK27538	ICLK27546
	3.0 mm	ICLK3008	ICLK3012	ICLK3016	ICLK3020	ICLK3025	ICLK3031	ICLK3038	ICLK3046
	3.5mm	ICLK3508	ICLK3512	ICLK3516	ICLK3520	ICLK3525	ICLK3531	ICLK3538	ICLK3546
	4.0mm	ICLK4008	ICLK4012	ICLK4016	ICLK4020	ICLK4025	ICLK4031	ICLK4038	-
	4.5mm	-	ICLK4512	ICLK4516	ICLK4520	ICLK4525	ICLK4531	-	-

TABLE 2		Balloon diameter (mm)							
		2.25	2.50	2.75	3.00	3.50	4.00	4.50	
Inflation pressure (atm)	6	2.08	2.31	2.54	2.77	3.26	3.74	4.15	
	7	2.13	2.37	2.61	2.85	3.35	3.84	4.27	
	8	2.19	2.43	2.68	2.93	3.43	3.92	4.38	
	9	2.25	2.50	2.75	3.00	3.50	4.00	4.50	
	10	2.31	2.55	2.83	3.08	3.58	4.09	4.61	
	11	2.36	2.61	2.90	3.15	3.66	4.17	4.68	
	12	2.40	2.65	2.95	3.20	3.71	4.25	4.74	
	13	2.44	2.70	3.01	3.25	3.77	4.32	4.80	
	14	2.48	2.74	3.06	3.29	3.82	4.37	4.86	
	15	2.51	2.77	3.09	3.32	3.86	4.41	4.92	
	16	2.55	2.80	3.12	3.36	3.90	4.46	4.98	
	17	2.59	2.83	3.15	3.39	3.94	4.50	5.04	
	18	2.62	2.86	3.19	3.43	3.98	4.54	5.11	
	19	2.65	2.89	3.23	3.47	4.03	4.59	5.19	
	20	2.68	2.93	3.27	3.51	4.08	4.65	5.27	
21	2.71	2.96	3.30	3.55	4.12	4.69	5.33		
22	2.74	2.99	3.34	3.60	4.18	4.74	5.40		

Pressure above Rated Burst Pressure
Do not exceed Rated Burst Pressure

These values result from in vitro testing

If postdilatation is required, a PTA balloon can be used. Use extreme care:
not to expand the 2.25 mm nominal diameter stents over 2.55 mm
not to expand the 2.50 mm and the 2.75 mm nominal diameter stents over 3.05 mm
not to expand the 3.00 mm and the 3.50 mm nominal diameter stents over 3.85 mm
not to expand the 4.00 mm and the 4.50 mm nominal diameter stents over 5.05 mm