
**CiD**  
Carbostent & Implantable Devices

INTERVENTIONAL CARDIOLOGY

CARBOFILM™ TECHNOLOGY FOR A BIO INDUCER SURFACE

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


**AVANTGARDE****„Carbofilm™ coated coronary stent on rapid exchange balloon catheter****1. DESCRIPTION**

The "AVANTGARDE CHRONO™ CARBOSTENT™" device consists of a coronary stent firmly held on the distal end of a semicompliant balloon catheter.

The **coronary stent** is a flexible implantable device that can be expanded using a PTCA catheter.

The stent is made of cobalt chromium alloy coated with „Carbofilm™, a thin turbostratic carbon film with a high-density crystalline structure substantially identical to that of the pyrolytic carbon used for mechanical cardiac valve discs. Coating the substrate with „Carbofilm™ provides it with the bio- and haemocompatible characteristics proper of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Two radio-opaque platinum markers are located at either end of the stent and enable accurate positioning of the device over the lesion to be treated.

The rapid-exchange type balloon catheter provides a safe means for conveying the coronary stent to the lesion to be treated.

The distal portion of the catheter has a hydrophilic coating and consists of two lumens: one is used for balloon inflation and deflation, the other for the passage of the guidewire.

Two radio-opaque markers, located outside the usable length of the balloon enable accurate placement across the stenosis.

The proximal portion of the catheter, with stainless steel hypotube, forms the lumen for balloon inflation and deflation.

Two depth indicators placed 90 and 100 cm from the distal end for femoral and brachial approach, respectively, aid in determining when the balloon is exiting the guiding catheter.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

The Manufacturer directly manufactures the AVANTGARDE CHRONO™ CARBOSTENT™ device and carries out all quality controls both during production and on the finished product, in accordance to the appropriate good manufacturing standards.

**How supplied:**

The AVANTGARDE CHRONO™ CARBOSTENT™ is supplied sterile, non-pyrogenic and individually packed in a pouch that **must not be placed in a sterile field**. The sterilisation made by the manufacturer uses an ethylene oxide and CO<sub>2</sub> mixture.

Sterility is guaranteed while the package remains intact correctly stored until the expiry date printed on the packaging (EXPIRY DATE).

**Contents:**

- One AVANTGARDE CHRONO™ CARBOSTENT™ device
- One compliance table

**Storage:**

Store in a cool, dry place, away from sunlight.

**2. INTENDED USE**

Use of the coronary stent is indicated in the treatment of coronary occlusions in order to maintain patency of the vessel.

**3. INDICATIONS**

The stent is indicated for use in the following cases:

- Elective treatment of obstructive lesions of native coronaries and aortocoronary bypasses in patients having the requisites for percutaneous transluminal coronary angioplasty (PTCA)
- When short- and long-term results of PTCA procedures are inadequate
- Primary angioplasty or intervention during acute myocardial infarction.

**4. CONTRAINDICATIONS**

The use of the stent is contraindicated in the following situations:

- Lesions considered untreatable with PTCA or other operating techniques.
- Disorders / allergies that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to contrast medium used during the procedure.
- Lesions on a vessel with a reference diameter < 2.25 mm.

**5. WARNINGS**

- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- **In patients with known hypersensitivity or allergies to metal components of the stent, the AVANTGARDE CHRONO™ CARBOSTENT™ stent should be used only if potential benefits of the stent outweigh potential risks.**
- There are no adequate or well-controlled studies in pregnant women for this product. The AVANTGARDE CHRONO™ CARBOSTENT™ stent should be used during pregnancy only if potential benefits of the stent outweigh potential risks.
- The device is intended for single use only. Do not reuse, reprocess or re-sterilize. This could compromise device performance and cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The AVANTGARDE CHRONO™ CARBOSTENT™ should be carefully handled

so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce alterations.

- The stent should not be touched. This is most important during catheter removal from packaging, placement over the guide wire, and advancement through the rotating hemostatic valve adapter and guiding catheter hub.
- The AVANTGARDE CHRONO™ CARBOSTENT™ device is designed to perform as a complete system. Do not use its components separately.
- The AVANTGARDE CHRONO™ CARBOSTENT™ device is also indicated for use in combination with PTCA. The delivery catheter is not a coronary dilation catheter; it should be used for stent deployment only.
- Do not use the AVANTGARDE CHRONO™ CARBOSTENT™ device if its packaging has been opened or damaged or if the stated "Expiry date" has passed. In such cases product performances and sterility are not guaranteed.
- Do not use the AVANTGARDE CHRONO™ CARBOSTENT™ device if the shaft of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- A heart surgery team should be available for a possible intervention.
- Follow the manufacturer's instructions for the use of accessories (guide catheter, guidewire, haemostatic valve).
- The AVANTGARDE CHRONO™ CARBOSTENT™ device should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- The traditional procedure for stent deployment requires pre-dilatation of the lesion.
- Recent scientific literature describes procedures performed without pre-dilatation.
- The stent should be implanted in the target lesion using its delivery catheter.
- If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the AVANTGARDE CHRONO™ CARBOSTENT™ device and the guide catheter as a single unit. Applying excessive force and/or incorrectly handling the system may lead to stent deployment or damage to the delivery catheter
- If the patient has multiple lesions in a single vessel, it is advisable to first treat the distal lesion, then the proximal lesions. This order will reduce the need to cross the proximal stent during implantation of the distal stent and therefore reduces the risk of displacing the proximal stent.
- Implantation of a coronary stent may cause dissection of the vessel distal and/or proximal to the stent and could also cause acute occlusion of the vessel, making a further operation necessary (CABG, further dilation, placement of additional stents or other procedures).
- Do not pressurise the system until the stent is positioned across the lesion to be treated
- Never use air or any other gas substances to inflate the balloon
- During catheter inflation, do not exceed the Rated Burst Pressure
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. This action may cause severe vessel damage.
- If the stent is lost in the coronary vessel, recovery procedures can be initiated. The procedures can, however, cause injury to coronary vessels and/or to the vessel access site.
- Do not attempt to clean or re-sterilize devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as dangerous medical waste with risk of infection

**6. PRECAUTIONS**

- Before the angioplasty procedure, examine the catheter to check that it operates correctly (no kinks, bends or other damage) and ensure that the dimensions are correct for the specific procedure in which it is to be used.
- Use of a gauged inflation device is strongly recommended.
- The AVANTGARDE CHRONO™ CARBOSTENT™ device should be used only by physicians trained specifically to perform percutaneous transluminal coronary angioplasty (PTCA) and coronary stent implantation.
- In cases of acute myocardial infarction, implantation of the stent should be carefully assessed by the operator due to the risk of acute thrombosis.
- Where vessel are excessively tortuous and there is also proximal atherosclerosis, the catheter may be difficult to advance. In such cases, incorrect handling could cause lead to dissection or rupture of the coronary vessel.
- Placement of a stent has the potential to compromise side branch patency.
- An unexpanded stent may be retracted into the guiding catheter one time only. An unexpanded stent should not be reintroduced into the artery once it has been pulled back into the guiding catheter.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.
- When multiple stents are needed, the stent materials must be of similar composition.

**7. MRI SAFETY INFORMATION**

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 2,890 G/cm (28.90 T/m)

**MRI Heating**

Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
<i>MR system reported, whole body averaged SAR</i>	2.9 W/kg	2.9 W/kg
<i>Calorimetry measured values, whole body averaged SAR</i>	2.1 W/kg	2.7 W/kg
<i>Highest temperature change</i>	2.4 °C	2.9 °C

**MR Artifact**

In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

## 8. DRUG SYSTEM

Studies in the clinical literature indicate the need to administrate appropriate anticoagulant anti-platelet and vasodilation therapy for successful stent implantation.

Antiplatelet treatment during the period after the procedure should be used per the guidelines from current ESC (2017 and subsequent updating) or ACC/AHA/SCAI Guidelines for PCI.

Given the complexity of the procedure and the many factors influencing the patient's condition, definition of the appropriate therapy is left to the discretion of the physician.

## 9. POTENTIAL ADVERSE EVENTS

Possible adverse effects of PTCA include, but are not limited to:

- Access site haematoma or pain
- Acute myocardial infarction
- Allergic reaction (to contrast medium, drugs used during the procedure or the materials composing the stent)
- Arterial spasm
- Arteriovenous Fistula
- Artery aneurism or pseudoaneurysm
- Cardiac arrhythmia
- Death
- Dissection, perforation, rupture of the coronary artery
- Distal embolism
- Embolization (air, tissue, device or thrombus)
- Haemorrhage
- Hypo/hypertension
- Infections
- Restenosis of the vessel
- Stent occlusion
- Thrombosis (acute, subacute or late)
- Unstable angina
- Ventricular fibrillation
- Vessel occlusion

## 10. DIRECTIONS FOR USE

### Auxiliary material (not included in the box)

- Guiding catheter(s) with a 5F diameter (inner diameter 1.47 mm) or larger
- Syringes
- Guidewires 0.014 inches in diameter (0.356 mm) or less
- Inflation device
- Introducer

### MODEL

Each "AVANTGARDE CHRONO™ CARBOSTENT" device is identified by a model code and a lot number; available models are listed in Table 1.

The code is composed of the letters AG, followed by two or three digits indicating the nominal diameter of expansion of the stent, and by two other digits indicating the stent length.

The lot number allows traceability of all information regarding the device manufacturing process and system control in the Quality Assurance archives of the Manufacturer.

To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included in each box; these labels can be attached to the patient's medical chart.

### Stent + delivery catheter preparation

The stent should be implanted in the target lesion using its delivery catheter.

- After having examined the package for any damage, remove the AVANTGARDE CHRONO™ CARBOSTENT device and bring it into a sterile area.
- Check that the delivery catheter does not present kinks, bends or other damage.
- Carefully slide the protective cover off the stent by grasping the cover at the distal end, do not touch the part of the cover over or proximal to the stent. Check that the stent is intact and well centered of the balloon.

**Damage to the AVANTGARDE CHRONO™ CARBOSTENT device may impair its performance. If the stent is out of place or damaged, do not use it.**

- Flush guidewire lumen with a heparin/saline mixture.
- Prepare the inflation device according to the manufacturer's instructions with dilute contrast medium.
- Purge any air from the balloon the stent is mounted, as follows:
  - 1) Fill the inflation device/syringe with 4 ml contrast medium, diluted (50% contrast medium and 50% sterile saline solution).
  - 2) After connecting the inflation device/syringe to the delivery catheter Luer connector, point the distal tip of the latter (balloon) downwards vertically.
  - 3) Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to neutral** as the system fills with contrast medium.
  - 4) Disconnect the inflation device/ syringe from the luer connection and remove all air.
  - 5) Without introducing air, repeat steps 2 and 3 and aspirate for 10-15 seconds until air bubbles no longer appear.

### Stent insertion

**The guiding catheters with a 5F (inner diameter 1.47 mm) or larger guide catheters are compatible for use with AVANTGARDE CHRONO™ CARBOSTENT devices.**

**Coronary guidewires with diameter of 0.014 inches (0.356 mm) or less**

**are compatible for use with the AVANTGARDE CHRONO™ CARBOSTENT device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.**

- Keep the delivery catheter the stent is crimped on **at ambient pressure during insertion** (neutral pressure).
- Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- Ensure that the haemostasis valve is completely open before inserting the system.
- Thread the distal tip of the catheter over the proximal end of the guidewire until the guidewire exits the catheter body at the guidewire exit port.
- Slowly advance the dilatation catheter until the appropriate depth marker aligns with the haemostatic valve hub. This indicates that the balloon is exiting the guiding catheter.
- Under fluoroscopy, slowly advance the system over the guidewire to allow the blood the retrograde filling of the guiding catheter; keep advancing the system until the stent reaches the treatment site.
- When advancing the delivery system into the target vessel, make sure that the stent and the delivery system move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon radiopaque markers in relation to the stent radiopaque markers.

**If resistance is encountered at any time, do not force the system: withdraw the AVANTGARDE CHRONO™ CARBOSTENT device and the guiding catheter as a single unit. Applying excessive force and/or incorrect handling may lead to stent deployment or damage to the delivery catheter.**

### Stent implant and expansion

- Position the stent and, with the help of the balloon's radiopaque markers, make sure it is correctly positioned in relation to the lesion to be treated and completely covering it.
- Expand the stent by slowly inflating the balloon at least to the nominal diameter. Table 2 lists the diameters of the balloon as inflation pressure varies (column I) for the seven nominal diameter groups: 2.25 mm (column II), 2.50 mm (column III), 2.75 mm (column IV), 3.00 mm (column V), 3.50 mm (column VI), 4.00 mm (column VII) and 4.50 mm (column VIII). The diameters shown on a dark background refer to pressures above the Rated Burst Pressure\*. The stent expands at a balloon inflation pressure of about 6 atm.
- Stent elastic recoil is between 2 and 7%\* depending on stent model and expansion diameter.**
  - \* NOTE: These values are reported from results of in-vitro testing.
- Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, reinflate the balloon for optimal implantation.**
- Before withdrawing the delivery catheter, maintain negative pressure for at least 30 seconds, as per routine PTCA procedure. Ensure the balloon is fully deflated before and remove the catheter. A device with larger and longer balloon requires a longer deflation time.
- If post-dilatation of the stent is required, a PTCA balloon can be used. Use extreme care:
  - not to expand the 2.25 mm nominal diameter stent over 2.55 mm;**
  - not to expand the 2.50 – 2.75 mm nominal diameter stents over 3.05 mm;**
  - not to expand the 3.00 -3.50 mm nominal diameter stents over 3.85 mm;**
  - not to expand the 4.00 – 4.50 mm nominal diameter stents over 5.05 mm;**

Optimal expansion requires the stent to be in full contact with the arterial wall, so that the stent diameter equals the diameter of the reference vessel.

**Ensure that the stent is fully expanded.**

## 11. LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packed with the utmost care, using the techniques considered most appropriate among those available in the current state of technology and applying integrated safety standards in the design and manufacture that will guarantee its safe use under the conditions and for the purposes stipulated, observing the precautions described in the paragraphs above, and that will nevertheless reduce, as far as possible, but not eliminate completely, the risks related to the use of the device.

The device must be used only under a specialist physician's responsibility, and taking into account any residual risks or possible side effects and complications of the treatment for which it is intended, among those described in other sections of this instruction booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the successful result following the use of the device or for its effectiveness in resolving a patient's condition. In fact, the results, in terms of both the patient's clinical status of the patient and the functionality and the working life of the device, depend on many factors beyond the manufacturer's control, among which are the patient's conditions, the surgical implant procedure or use, and the methods of handling the device after the pack has been opened.

In the light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. To this end, the purchaser must return the device to the Manufacturer, which reserves the right to inspect the device deemed to be faulty and to determine whether the device actually has manufacturing or material defects. The warranty consists exclusively of the replacement of the device found to be faulty with another device of the same or similar type.

The warranty only applies if the device is returned correctly packed to the Manufacturer, and is accompanied by a detailed report in writing describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the purchaser the expenses incurred for replacement of the device acknowledged to be faulty. Manufacturer does not assume any liability in any case of failure to observe the

methods of use and the precautions stated in this instruction booklet and in the event of use of the device after the expiry date stated on the pack. Furthermore, the Manufacturer does not assume any liability in relation to the consequences arising from the medical choices and the methods of use or the application of the device. The Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, as a result of the application of the device, and the choice of implant technique used by the operator. The agents and the representatives of the Manufacturer are not authorized to amend any of the conditions of this warranty or to undertake any further commitments or to offer any guarantee in relation to this product, beyond the terms stated above.

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
TABLE 1		STENT LENGTH (mm)							
		7mm	8mm	12mm	16mm	20mm	24mm	25mm	31mm
NOMINAL STENT DIAMETER (mm)	2.25 mm	ICV9601	-	ICV9602	ICV9603	ICV9604	ICV9605	-	-
	2.5 mm	-	ICV9606	ICV9607	ICV9608	ICV9609	-	ICV9610	-
	2.75 mm	-	ICV9611	ICV9612	ICV9613	ICV9614	-	ICV9615	-
	3.0 mm	-	ICV9616	ICV9617	ICV9618	ICV9619	-	ICV9620	ICV9621
	3.5 mm	-	ICV9622	ICV9623	ICV9624	ICV9625	-	ICV9626	ICV9627
	4.0 mm	-	ICV9628	ICV9629	ICV9630	ICV9631	-	ICV9632	ICV9633
	4.5 mm	-	-	ICV9634	ICV9635	ICV9636	-	ICV9637	ICV9638

TABLE 2		BALLOON DIAMETER (mm)						
Inflation pressure (atm)								
	2.25	2.50	2.75	3.00	3.50	4.00	4.50	
6	2.08	2.31	2.54	2.77	3.26	3.74	4.15	
7	2.13	2.37	2.61	2.85	3.35	3.84	4.27	
8	2.19	2.43	2.68	2.93	3.43	3.92	4.38	
<b>9</b>	<b>2.25</b>	<b>2.50</b>	<b>2.75</b>	<b>3.00</b>	<b>3.50</b>	<b>4.00</b>	<b>4.50</b>	
10	2.31	2.55	2.83	3.08	3.58	4.09	4.61	
11	2.36	2.61	2.90	3.15	3.66	4.17	4.68	
12	2.40	2.65	2.95	3.20	3.71	4.25	4.74	
13	2.44	2.70	3.01	3.25	3.77	4.32	4.80	
14	2.48	2.74	3.06	3.29	3.82	4.37	4.86	
15	2.51	2.77	3.09	3.32	3.86	4.41	4.92	
16	2.55	2.80	3.12	3.36	3.90	4.46	4.98	
17	2.59	2.83	3.15	3.39	3.94	4.50	5.04	
18	2.62	2.86	3.19	3.43	3.98	4.54	5.11	
19	2.65	2.89	3.23	3.47	4.03	4.59	5.19	
20	2.68	2.93	3.27	3.51	4.08	4.65	5.27	
21	2.71	2.96	3.30	3.55	4.12	4.69	5.33	
22	2.74	2.99	3.34	3.60	4.18	4.74	5.40	

Pressure above Rated Burst Pressure  
Do not exceed Rated Burst Pressure

*These values result from in vitro testing*

If postdilatation is required, a PTCA balloon can be used. Use extreme care:  
**not to expand the 2.25 mm nominal diameter stent over 2.55 mm**  
**not to expand the 2.50 mm and the 2.75 mm nominal diameter stents over 3.05 mm**  
**not to expand the 3.00 mm and the 3.50 mm nominal diameter stents over 3.85 mm**  
**not to expand the 4.00 mm and the 4.50 mm nominal diameter stents over 5.05 mm**