INTERVENTIONAL CARDIOLOGY

CRE8 EVO™

AMPHILIMUS™ (SIROLIMUS+FATTY ACID)
ELUTING CORONARY STENT






abLUMINAL
réservoir
technology 

IC0730.137.102 issued 201910 Rev.4

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of  **ALVIMEDICA** group
www.alvimedica.com

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CRE8™ EVO
Sirolimus eluting μ Carbofilm™ coated coronary stent on rapid exchange balloon catheter

1. DESCRIPTION

The CRE8™ EVO device consists of a Sirolimus eluting coronary stent firmly held on the distal end of a semicompliant balloon catheter.

The coronary **stent** is a flexible implantable device that can be expanded using a PTCA catheter.

The stent is made of cobalt chromium alloy (L605) and is coated with μ Carbofilm™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs. Coating the substrate with μ Carbofilm™ provides it with the bio- and haemocompatible characteristics proper to pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

The outer surface of the stent has dedicated grooves, fully coated with μ Carbofilm™, for containing the pharmaceutical formulation Amphilimus™, which is composed of the drug Sirolimus and a mixture of long-chain fatty acids.

The specific drug dosage for the treated vessel surface is approximately 0.9 $\mu\text{g}/\text{mm}^2$, which corresponds to a minimum dose of 50 μg on the smaller stent (2.0x9 mm) and a maximum dose of 395 μg on the larger stent (3.5x46 mm).

Two radio-opaque platinum markers at either end of the stent enable accurate positioning on the lesion to be treated.

The rapid-exchange type **balloon catheter** provides a safe means for delivering the coronary stent to the lesion to be treated.

The distal portion of the catheter consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radio-opaque markers, located outside the usable length of the balloon, enable accurate placement across the stenosis.

The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.

Two indicators placed at a depth of 90 and 100 cm from the distal end, signal the end of the balloon catheter from the guiding catheter, in cases of brachial or femoral approach respectively.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

The Manufacturer directly manufactures the CRE8™ EVO device and carries out all quality controls both during production and on the finished product, according to the norms of Good Manufacturing Practice.

2. INTENDED USE

The coronary stent is indicated for the treatment of stenotic lesions in coronary arteries and coronary artery bypasses to maintain vessel patency.

Randomized clinical trials have shown that drug-eluting stents can significantly reduce angiographic late loss (late lumen loss), binary restenosis and repetition of target lesion revascularization.

3. INDICATIONS

The stent is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease due to de novo or restenotic lesions of native coronary arteries even in presence of:

- Acute coronary syndrome,
- Concomitant diabetes
- Multiple lesions

4. CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations.

- Lesions considered untreatable with PTCA or other interventional techniques
- Patients who do not qualify for coronary artery bypass surgery
- Women who are pregnant
- Disorders/allergies that limit the use of anti-platelet and/or anticoagulant therapy
- Severe allergy to the contrast medium or drugs used during the procedure
- Stenosis of an unprotected vessel
- Ejection fractions $\leq 30\%$
- Diffuse distal coronary artery disease
- Lesions in an area involving a major side branch. A major side branch is defined as a vessel that would be bypassed if it was occluded.
- Lesions on a vessel with reference diameter < 2.0 mm
- Lesions distal to stenosis $\geq 50\%$ non dilatable
- Patients with known hypersensitivity or allergies to Sirolimus, fatty acids (such as stearic acid, palmitic acid, behenic acid) or the metal components of the stent

5. MODEL

Each CRE8™ EVO device is identified by a model code and a lot number; available product codes are listed in Table 1. The code is composed of the letters IC, followed by two letters identifying the type of device, two or three digits indicating the nominal diameter of expansion of the stent and a further two digits indicating the length of the stent. The batch number allows traceability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included on each box; these labels can be attached to the implant patient's medical chart.

6. HOW SUPPLIED

The CRE8™ EVO device is supplied sterile in a foil pouch that must not be brought into a sterile area.

NOTE: The foil pouch is the only sterile barrier.

The manufacture sterilisation uses an ethylene oxide and CO₂ mixture.

Sterility and chemical stability are guaranteed as long as the package remains intact and correctly stored, until the expiry date printed on the packaging (USE BEFORE DATE).

7. STORAGE

Store in a cool, dry place at a temperature of 25°C; excursions permitted to 30°C.

8. WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause risks of device contamination and infection for the patient, inflammation and patient-to-patient transmission of infective diseases.
- The CRE8™ EVO device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce alterations.
- The stent should not be touched with bare hands.
- **The stent should not come into contact with liquids prior to preparation and positioning. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).**
- The CRE8™ EVO device is designed to perform as a system. Do not use its components separately.
- The CRE8™ EVO device is indicated for use in combination with PTCA. The delivery catheter is not a coronary dilatation catheter; it should be used for stent deployment only.
- Do not use the CRE8™ EVO device if its packaging has been opened or damaged, incorrectly stored, or if it is past the stated expiry date. In such cases product sterility is not guaranteed.
- Do not use the CRE8™ EVO if, during handling, due to torsion or resistance to insertion, the proximal portion of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The CRE8™ EVO device should be used only by physicians trained specifically to perform percutaneous transluminal coronary angioplasty (PTCA) and coronary stent implantation.
- A heart surgery team must be available for a possible intervention.
- The CRE8™ EVO device should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- The traditional procedure for stent deployment requires pre-dilatation of the lesion.
- Recent scientific literature describes procedures performed without pre-dilatation.
- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- The stent should be implanted in the target lesion using its delivery catheter.
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
- Do not attempt to clean or re-sterilise devices that have been in contact with blood and organic tissues. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- Follow the manufacturer's instructions when using accessories (guiding catheter, guidewire, haemostatic valve).
- If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the CRE8™ EVO device and the guiding catheter as a single unit. The application of excessive force and/or incorrect handling of the system may lead to stent deployment or damage to the delivery catheter.
- If the patient has multiple lesions in a single vessel, it is advisable to first treat the distal lesion, then the proximal lesions. This order will reduce the need to cross the proximal stent during implantation of the distal stent and therefore reduces the risk of displacing the proximal stent.
- Implantation of a coronary stent may cause dissection of the vessel distal and/or proximal to implantation site and may also cause acute occlusion of the vessel, making an additional intervention necessary (CABG, further dilatation, placement of additional stents, or other procedures).
- If the stent is lost in the coronary vessel, recovery procedures can be initiated. The procedures can, however, cause injury to coronary vessels and/or to the vessel access site.
- During catheter inflation do not exceed the Rated Burst Pressure.
- Never use air or any other gas to inflate the balloon.

9. PRECAUTIONS

- Use of a gauged inflation device is strongly recommended.
- In cases of acute myocardial infarction, implantation of the stent should be carefully assessed by the operator due to the risk of acute thrombosis.
- Where vessels are excessively tortuous and there is also proximal atherosclerosis, the catheter may be difficult to advance. In such cases incorrect handling could lead to dissection or rupture of the coronary vessel.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.

10. MRI SAFETY INFORMATION

Non-clinical testing demonstrated that the stent is MR Conditional.

A patient with this device can be scanned safely, immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
- Highest spatial gradient magnetic field of 720-Gauss/cm or less
- Maximum MR system reported, whole body averaged SAR of 2.9-W/kg for 15-min of imaging

In non-clinical testing, the stent under MRI scanning for 15-min in 1.5-Tesla/64-MHz is expected to produce a maximum temperature increase of 2.4°C and under MRI scanning for 15-min in 3-Tesla/128-MHz a maximum temperature increase of 2.9°C. The temperature changes will not pose a hazard to patients under the conditions indicated. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent.

11. DRUG REGIME

Studies in clinical literature indicate the need to administer anticoagulant therapy during the procedure and anti-platelet treatment after the procedure.

Antiplatelet treatment during the period after the procedure involves the administration of:

- Thienopyridine (ticlopidine or clopidogrel or prasugrel in the doses indicated in the relative instruction leaflets or according to medical prescription)
- Aspirin at a minimum amount of 100 mg/day indefinitely.

12. INTERACTIONS WITH DRUGS

Although no specific clinical data are available, certain drugs, such as Tacrolimus, which act through the same binding protein (FKBP), may interfere with the efficacy of Sirolimus. There have been no studies on the interaction between drugs. Sirolimus is metabolized by CYP3A4. Strong CYP3A4 inhibitors may cause greater exposure to Sirolimus up to levels associated with systemic effects, especially in the case of multiple implants. Systemic exposure to Sirolimus if the patient is taking concomitant systemic immunosuppressive drugs must also be considered.

13. DIRECTIONS FOR USE

Auxiliary material (not included in the package)

- Guiding catheter(s) with a 5F diameter (inner diameter 1.47 mm) or larger
- Syringes
- Guidewires 0.014 inches in diameter (0.356 mm) or less
- Inflation device
- Introducer.

Stent + delivery catheter preparation

Before commencing the angioplasty procedure, administer heparin and check that the patient's Activated Clotting Time (ACT) is greater than 300 seconds.

During preparation of the system, take special care not to let to any liquid come into contact with the stent. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline solution, the contact time should be limited (maximum one minute).

The stent should be implanted in the target lesion using its delivery catheter. After having examined the package for any damage, remove the CRE8™ EVO device and bring it into a sterile area.

- Check that the delivery catheter does not present kinks, bends or other damage.
- Carefully slide the protective cover off the stent by grasping the cover at the distal end. Check that the stent is intact and well centred on the balloon.

Damage to the CRE8™ EVO device may impair its performance.

- Flush guidewire lumen with heparin/saline mixture. **Caution: do not place the stent in contact with liquids while flushing the guidewire lumen.**
- Prepare the inflation device according to the manufacturer's instructions.
- Purge any air from the balloon the stent is mounted on as follows:
 - Fill the inflation device with 4 ml of contrast medium;
 - After connecting the inflation device to the delivery catheter Luer connector, point the distal tip of the catheter (balloon) downwards vertically;
 - Apply negative pressure and aspirate for at least 30 seconds. **Allow the pressure to rise gradually back to normal** as the system fills with contrast medium.
 - Without introducing air, repeat step 3 and aspirate for 10-15 seconds until air bubbles no longer appear.

Stent insertion

The guiding catheters with a 5F diameter (inner diameter 1.47 mm) or larger are compatible for use with the CRE8™ EVO device.

Coronary guidewires 0.014 inches in diameter (0.356 mm) or less are compatible for use with the CRE8™ EVO device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- Keep the delivery catheter the stent is crimped on at ambient pressure during insertion.
- Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- Ensure that the haemostasis valve is completely open before introducing the system.
- Slowly advance the system over the guidewire to allow the blood the retrograde filling of the guiding catheter; keep advancing the system until the stent reaches the treatment site.
- When advancing the delivery system into the target vessel, make sure that the stent and the delivery balloon move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon's radiopaque markers in relation to the stent's radiopaque markers.

If resistance is encountered at any time, do not force the system: withdraw the CRE8™ EVO device and the guiding catheter as a single unit. The application of excessive force and/or incorrect handling of the system may lead to stent deployment or damage to the delivery catheter.

Stent implantation and expansion

- Position the stent and, with the help of the **balloon's radiopaque markers**, make sure it is correctly positioned in relation to the lesion to be treated and completely covered.
- Expand the stent by slowly inflating the balloon at least to the nominal diameter. Table 2 lists the diameters of the balloon as inflation pressure varies (column I), for the eight nominal diameter groups: 2.00 mm (column II), 2.25 mm (column III), 2.5 mm (column IV), 2.75 mm (column V), 3.0 mm (column VI), 3.5 mm (column VII), 4.0 mm (column VIII) and 4.5 mm (column IX). Values on dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a minimum balloon inflation pressure of about 5 atm.
- Stent elastic recoil is less than 7%* depending on stent model and expansion diameter.**
* NOTE: These values are obtained from the results of in-vitro testing.
- Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation.**

- Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTCA procedure.
- If post-dilatation of the stent is required, a PTCA balloon can be used. Use extreme care:

**not to expand the 2.00 mm nominal diameter stent over 3.05 mm.
not to expand the 2.25 mm nominal diameter stent over 3.05 mm.
not to expand the 2.50 mm nominal diameter stent over 3.85 mm.
not to expand the 2.75 mm nominal diameter stent over 3.85 mm.
not to expand the 3.00 mm nominal diameter stent over 3.85 mm.
not to expand the 3.50 mm nominal diameter stent over 5.05 mm.
not to expand the 4.00 mm nominal diameter stent over 5.05 mm.
not to expand the 4.50 mm nominal diameter stent over 5.05 mm.**

Optimal expansion requires the stent to be in full contact with the arterial wall, so that the stent diameter equals the diameter of the reference vessel. ENSURE THAT THE STENT IS FULLY EXPANDED.

14. ADVERSE EVENTS

Coronary stent implantation may give rise to the following complications:

- Access site haematoma
- Acute myocardial infarction
- Allergic reaction (to contrast medium, drugs used during the procedure or the materials composing the stent)
- Arterial spasm
- Artery aneurysm or pseudoaneurysm
- Cardiac arrhythmia
- Death
- Dissection, perforation, rupture of the artery
- Distal embolism
- Haemorrhage
- Infection and pain at the access site
- Restenosis of the vessel
- Stent embolization
- Stent occlusion
- Thrombosis (acute, subacute or late)
- Unstable angina
- Ventricular fibrillation
- Vessel occlusion

15. LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the Manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. Under such circumstances, the customer shall deliver the device to the Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty consists exclusively in the replacement of the defective device with another device of the same type or equivalent, from the same Manufacturer.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for return of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use-before date printed on the package.

Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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Pressure above Rated Burst Pressure
Do not exceed Rated Burst Pressure

These values result from in vitro testing

If postdilatation is required, a PTCA balloon can be used. Use extreme care:
not to expand the 2.00 mm and 2.25 mm nominal diameter stents over 3.05 mm
not to expand the 2.50 mm, 2.75 mm and 3.00 mm nominal diameter stents over 3.85 mm
not to expand the 3.50 mm, 4.00 mm and 4.50 mm nominal diameter stents over 5.05 mm