PTCA Rapid exchange balloon catheter

PTCA Rapid exchange balloon catheter
Catetere a scambio rapido per PTCA
Schnellwechsel-Katheter für PTCA
Cathéter à échange rapide pour ACTP
Catéter de balón de cambio rápido para ACTP
Kateter til PTCA til hurtig udskiftning
Cateter de balão de troca rápida para ACPT
Snel verwisselbare PTCA-ballonkatheter
Snabbt utbytbar ballongkateter för PTCA
Pikavaihdettava PCTA-pallokatetri
Καθετήρας ταχείας εναλλαγής για διαδερμική στεφανιαία αγγειορλαστική (PTCA)
PTCA hızlı değişşim balon kateteri
Rychle měnitelný balónkový katetr pro PTCA

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I00704-137-104 ISSUED 201908 Rev. 8

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of **Alvimedica** group
www.alvimedica.com

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DESCRIPTION

The PRO-HP PEGASO™ percutaneous transluminal coronary angioplasty (PTCA) catheter (shown in the figure 1) is a rapid exchange type catheter.

It has a low profile balloon (less than 1 mm in diameter before use). Balloon diameters at different pressures (compliance) are shown in the table 2.

The distal portion of the catheter consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radiopaque markers located at both ends of the balloon permit accurate placement of the balloon across the stenosis under fluoroscopy.

The proximal portion, a stainlees steel hypotube, contains the lumen for balloon inflation and deflation.

Two depth markers, one at 92 and the other at 102 cm from the distal tip, aid in determining when the balloon is exiting the guiding catheter.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

INDICATIONS

The PRO-HP PEGASO™ catheter is indicated for:

- patients who are candidates for coronary artery bypass graft (CABG) surgery who may benefit from a less invasive technique like PTCA;
- patients who have previously undergone CABG surgery and present with recurrent symptoms and progression of disease in the native coronary artery or stenosis and closure of the graft;
- patients suffering from atherosclerotic coronary artery disease with a short lesion that is concentric, discrete, and can be accessed by a dilatation catheter;
- patients with native atherosclerotic coronary disease. (The physician must determine the accessibility to each lesion separately and decide whether PTCA constitutes the therapy of choice.)

CONTRAINDICATIONS

The PRO-HP PEGASO™ catheter is contraindicated in patients presenting with one or more of the following problems:

- left main coronary artery disease;
- previously diagnosed coronary artery spasm in diseased coronary segments;
- ineligibility for coronary artery bypass surgery.

PACKAGING

The catheter is supplied sterile and individually packed in a pouch that must not be placed in a sterile field.

The sterilization made by the manufacturer uses an ethylene oxide and CO₂ mixture.

Sterility is guaranteed while the package remains intact and until the expiry date printed on the packaging (EXPIRY DATE).

PRECAUTIONS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The PTCA catheter should only be used by or under the direction of physicians thoroughly trained in the technique of catheter dilation or PTCA.
- The PTCA catheter should only be used when a full cardiac surgery team and operating room are on stand-by.
- The catheter must be handled under fluoroscopic control using X-ray equipment that produces high-quality images.
- Do not use the catheter if during handling, due to torsion or resistance to insertion, the proximal portion of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The PTCA catheter should only be advanced over a guidewire. Do not remove the guidewire from the catheter during use. Do not advance the catheter unless the balloon is completely deflated and under vacuum.
- The inflated balloon diameter should never exceed the diameter

of the coronary vessel in the area just proximal or distal to the stenosis.

- Do not attempt to reclean or resterilise devices that have been in contact with blood or body tissues. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- During catheter inflation do not exceed the Rated Burst Pressure, as excessive pressure may cause the balloon or the inflation device to malfunction. Use of a gauged inflation device is strongly recommended.
- The short-term and long-term effect at pressures above the nominal diameter inflation pressure are not known. Never use air or other gases to inflate the balloon.

COMPLICATIONS

Possible complications arising after PTCA include: prolonged angina, coronary occlusion, coronary spasm, dissection/intimal tear, hypotension, bradycardia, ventricular fibrillation, coronary embolism, myocardial infarction.

These complications could require emergency CABG surgery or other medical treatment, which could result in death.

INSTRUCTIONS FOR USE

Catheter preparation

Using sterile techniques, remove the catheter from its package without bending it, and bring it into a sterile field.

Examine the catheter for any damage.

Purge any air from the catheter as follows:

1. Fill a 20-30 cc syringe with 3-5 ml contrast mixture (50% contrast medium and 50% sterile saline solution and purge the remaining air.
2. Attach the syringe to the Luer connection.
3. With the tip of the catheter pointing downwards, aspirate for 30 seconds. Slowly release the plunger of the syringe and allow the contrast medium to flow into the distal portion of the catheter.
4. Repeat step 3, aspirate the device for 10-15 seconds until bubbles no longer appear. Remove the syringe.
5. Prepare the inflation device according to the manufacturer's instructions. Attach the inflation device to the inflation port to make fluid-to-fluid contact.
6. Before inserting the catheter, apply negative pressure to the balloon with the inflation device. The negative pressure will provide the lowest possible balloon profile to ease catheter insertion.
7. Carefully remove the distal stylet and the balloon cover.

Catheter insertion

1. In the standard manner prepare the catheter insertion site for percutaneous introduction. Cannulate the artery for catheter introduction.

2. Under fluoroscopic control:

- a) Place a guiding catheter in the aortic arch.

NOTE: All 5F (1.47 mm) or larger guiding catheters are compatible for use with PRO-HP PEGASO™ catheters.

- b) Using the bare-wire technique, place a guidewire across the stenosis to be treated.

NOTE: All guidewires 0.014 inches (0.356 mm) or less in diameter are compatible for use with PRO-HP PEGASO™ catheters. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

3. Thread the distal tip of the catheter over the proximal end of the guidewire until the guidewire exits the catheter body at the guidewire exit port.

4. Advance the catheter over the guidewire until the catheter radiopaque markers show the balloon is centered in the area to be dilated.

5. With the catheter positioned within the stenosis, inject the diluted contrast medium into the inflation lumen of the catheter to inflate the balloon.

6. The success of the dilatation procedure can be determined by injecting contrast medium through the guiding catheter.

NOTE: It is usually not necessary to try to restore a stenotic area to normal lumen diameter.

7. After dilatation, completely deflate the balloon and remove the dilatation catheter.

LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. The buyer shall return the device to the Manufacturer, which reserves the right to inspect the returned device and, at its sole discretion, to determine whether the device is truly defective in manufacturing or material. The warranty consists exclusively in the replacement of the defective device with a Manufacturer device of the same fabrication or equivalent.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for replacement of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package.

Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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