CiD
Carbostent & Implantable Devices

PERIPHERAL VASCULAR

CARBOSTENT™

RADIX2

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
IC0711.137.103 issued 202011 Rev. 4

CiD
Carbostent & Implantable Devices

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of **Alvimedica** group
www.alvimedica.com

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DESCRIPTION

The "CARBOSTENT™ RADIX2" renal stent system device consists of a renal stent firmly held on the distal end of a semicompliant balloon catheter.

The **peripheral renal stent** is a flexible implantable device that can be expanded using a PTRA balloon catheter.

The stent is made of AISI 316 LVM stainless steel and is coated with Carbofilm™, a thin carbon film with a high-density turbostratic structure substantially identical with that of the pyrolytic carbon used for mechanical cardiac valve discs.

Coating the substrate provides it with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Two radio-opaque platinum markers at either end of the stent enable accurate positioning over the lesion to be treated.

The rapid-exchange type **balloon catheter** provides a safe means for delivering the renal stent to the lesion to be treated.

The distal portion of the catheter consists of two lumens: one for balloon inflation and deflation, the other for guidewire advancement and retraction.

Two radio-opaque markers, located outside the usable length of the balloon, enable accurate placement across the stenosis.

The proximal portion of the catheter, a stainless steel hypotube, contains the lumen for balloon inflation and deflation.

The proximal end of the catheter has a female Luer lock for attachment to an inflation device.

The Manufacturer directly manufactures the "CARBOSTENT™ RADIX2" renal stent system and carries out all quality controls both during production and on the finished product, according to the norms of Good Manufacturing Practice.

INTENDED USE

The peripheral renal stent is intended for use in the treatment of atherosclerotic lesions of the renal arteries in order to maintain vessel patency.

INDICATIONS

The stent is specifically indicated for the treatment of occlusive lesions of the native renal arteries. The stent is also indicated in the case of inadequate results (residual stenosis > 30%) or dissection of the vessel wall following PTRA procedure with or without stenting.

CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations.

- Renal vessel to be treated with evidence of massive thrombosis prior to PTRA procedure.
- PTRA treatment technically unsuccessful, for example due to impossibility of accessing the site with the guiding catheter or inflation catheter.
- Aneurysm of the abdominal aorta involving the renal ostium and immediately adjacent to the stent implant site.
- Aneurysm of the renal arteries.
- Disorders that limit the use of anti-platelet and/or anticoagulant therapy.
- Severe allergy to pharmaceuticals and/or contrast medium commonly used in PTRA procedures.
- Lesions considered untreatable with PTRA or other interventional techniques.
- Lesions with severe calcification that resists dilatation with risk of vessel wall rupture.

MODEL

Each "CARBOSTENT™ RADIX2" renal stent system device is identified by a model code and a lot number; available models are listed in Table 1.

The model code consists of the letters RW, which specify the name of the renal stent system, two digits indicating the nominal expanded stent diameter and two digits indicating the stent length, followed by letters L or S according to whether the delivery catheter is 150 or 75 cm long.

The lot number allows trackability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate trackability of the device at the user end, the product code is printed on the adhesive labels included in each box; these labels can be attached to the implant patient's medical chart.

PACKAGING

The "CARBOSTENT™ RADIX2" renal stent system device is supplied sterile and individually packed in a pouch that **must not be brought into a sterile area**. The manufacture sterilisation uses an ethylene oxide and CO₂ mixture.

Sterility is guaranteed while the package remains intact and until the expiry date printed on the packaging (USE BEFORE DATE).

WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The "CARBOSTENT™ RADIX2" renal stent system device is designed to perform as a system. Do not use its components separately.
- The "CARBOSTENT™ RADIX2" renal stent system delivery catheter is not a peripheral dilatation catheter; it should be used for stent deployment only.
- Do not use the "CARBOSTENT™ RADIX2" renal stent system if its packaging has been opened or damaged or if the period of sterility has expired. In such cases product sterility is not guaranteed.
- The "CARBOSTENT™ RADIX2" renal stent system should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce mechanical alterations.

- Do not use the "CARBOSTENT™ RADIX2" renal stent system if, during handling, due to torsion or resistance to insertion, the proximal portion of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The "CARBOSTENT™ RADIX2" renal stent system should be used by physicians thoroughly trained in percutaneous transluminal renal angioplasty (PTRA) and in peripheral and/or coronary stent placement.
- The "CARBOSTENT™ RADIX2" renal stent system should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- Always advance the delivery catheter balloon with a guidewire.
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially deployed stent. Attempted repositioning may cause severe vessel damage.
- Do not attempt to clean or re-sterilise devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- Follow manufacturer's indications for the use of accessories (sheath and/or guiding catheter, guidewire, haemostatic valve).
- If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the system and the guiding catheter as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.
- Implantation of a renal stent may cause dissection of the vessel distal and/or proximal to implantation site and may also cause sudden occlusion of the vessel, making an additional intervention necessary (surgical vascular replacement operation, further dilatation, placement of additional stents, or other procedures).
- Extreme caution should be exercised to reduce the risk of plaque embolism.
- If the stent is lost in the vessel, recovery procedures can be initiated. The procedures can, however, cause injury to vessels and/or to the vessel access site.
- During catheter inflation do not exceed the Rated Burst Pressure.
- Never use air or any other gas to inflate the balloon.

PRECAUTIONS

- Use of a gauged inflation device is strongly recommended.
- When an unfavourable anatomy (vessel angulation) associated with proximal atherosclerosis is present, catheter advancement may be difficult. In such cases incorrect handling could cause dissection, perforation or rupture of the vessel.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.
- To minimise the risk of stent migration in the presence of a strong magnetic field, patients who have received a stent should not undergo magnetic resonance imaging (MRI) until the stent is fully endothelialised (ca. 8 weeks). The stent might affect the magnetic field, creating artefacts in the MRI images.

DRUG REGIME

Studies in clinical literature indicate the need to administer pre-operation anti-platelet therapy, anticoagulant therapy during the procedure and post-operation anti-platelet treatment.

Given the complexity of the procedure and the many factors influencing patient condition, definition of the appropriate therapy is left to the discretion of the physician.

DIRECTIONS FOR USE

Stent + delivery catheter preparation

Prior to stent delivery, administer Heparin and check that the patient's Activated Clotting Time (ACT) is greater than at least 180 sec.

The stent should be implanted in the target lesion using its delivery system.

- a) After having examined the package for any damage, remove the "CARBOSTENT™ RADIX2" renal stent system and bring it into a sterile area.
- b) Check that the delivery catheter does not present kinks, bends or other damage.
- c) Pull back the distal stylet.
- d) Carefully slide the protective cover off the stent by grasping the cover at the distal end, then the stylet. Check that the stent is intact and well centred on the balloon.

Damage to the "CARBOSTENT™ RADIX2" renal stent system may impair its performance.

- e) Flush balloon catheter guidewire lumen with heparin/saline mixture.
- f) Prepare the inflation device according to the manufacturer's instructions.
- g) Purge any air from the balloon the stent is mounted on as follows:
 - 1) Fill the inflation device with 4 ml of contrast mixture;
 - 2) After connecting the inflation device to the delivery catheter Luer connector, point the distal tip of the catheter downwards;
 - 3) Apply negative pressure and aspirate for at least 30 seconds. **Allow the**

pressure to rise gradually back to normal as the system fills with contrast medium.

- 4) Without introducing air, repeat step 3 and aspirate for 10-15 seconds until bubbles no longer appear.

Stent insertion

The "CARBOSTENT™ RADIX2" renal stent system device of diameter 5.0 and 6.0 mm is compatible with guiding catheters of diameter 6F (internal diameter 1.73 mm) or greater, and with sheaths of diameter 5F (internal diameter 1.67 mm) or greater.

The "CARBOSTENT™ RADIX2" renal stent system device of diameter 7.0 mm is compatible with 6F guiding catheters with large lumen (internal diameter 1.8 mm) or greater, and with 5F sheaths with large lumen (internal diameter 1.78 mm) or greater.

Guidewires 0.018 inches (0.47 mm) or less are compatible for use with the "CARBOSTENT™ RADIX2" renal stent system. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- a) Keep the delivery catheter the stent is pre-mounted on at ambient pressure during insertion.
- b) Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast medium.
- c) Ensure that the haemostasis valve is completely open before introducing the stent + delivery catheter system.
- d) Slowly advance the system over the guidewire until the stent reaches the treatment site.
- e) When advancing the delivery system into the target vessel, make sure that the stent and the delivery balloon move as a single unit. This should be done by observing, under fluoroscopy, the position of the balloon markers in relation to the stent markers.

If resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the system and the guiding catheter (and/or sheath) as a single unit. Applying excessive force and/or incorrectly handling the system may lead to loss of the stent or damage to the delivery catheter.

Stent deployment and expansion

- a) Position the stent across the target lesion.
- b) Deploy the stent by slowly inflating the balloon to the nominal diameter. Table 2 (pressure/diameter) lists the diameters of the system as balloon inflation pressure varies (column I), for the three nominal diameter groups: 5.0 mm (column II), 6.0 mm (column III) and 7.0 mm (column IV). Values on dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a balloon inflation pressure of about 6 atm.

Stent elastic recoil is between 2 and 5%* depending on stent model and expansion diameter.

*** NOTE: These values are obtained from the results of in-vitro testing**

- c) Deflate the balloon and check angiographically that the stent has fully expanded. If necessary, re-inflate the balloon for optimal implantation.
- d) Before withdrawing the delivery catheter, maintain negative pressure for a minimum of 30 seconds, as per routine PTRA procedure.
- e) If post-dilatation is required, a PTRA balloon can be used. Use extreme care:
not to expand the 5.0 mm nominal diameter stent over 6.5 mm
not to expand the 6.0 mm nominal diameter stent over 6.5 mm
not to expand the 7.0 mm nominal diameter stent over 8.0 mm.

The final diameter of the stent must be suitable for the reference vessel diameter.

Ensure that the stent is fully expanded.

COMPLICATIONS

Stent implantation may give rise to the following complications:

- Death
- Renal necrosis
- Dissection, perforation, rupture of the renal artery and abdominal aorta
- Perforation of the renal parenchyma
- Stent embolism
- Chronic stent occlusion (restenosis)
- Vessel occlusion
- Arterial spasm
- Infection
- Access site haematoma
- Drug and contrast medium induced collateral effects (anti-coagulant/anti-platelet)
- Acute or sub-acute thrombosis localised at stent level

LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforesaid conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its intended use, as mentioned in other sections of this instructions booklet.

Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of

the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. The buyer shall return the device to the Manufacturer, which reserves the right to inspect the returned device and, at its sole discretion, to determine whether the device is truly defective in manufacturing or material. The warranty consists exclusively in the replacement of the defective device with a device of the same fabrication or equivalent.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for replacement of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package. Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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