
EASY FLYPE & HIFLYPE

CARBOSTENT™

100721.137.101 Issued 202105 Rev.7

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of  Alvimedica group
www.alvimedica.com

EASY FLYPE & HIFLYPE

CARBOSTENT™

BIO INDUCER SURFACE
SELF EXPANDING STENT

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



Figure 1

Detail B

1	Radiopaque tip	9	Thumbwheel
2	Guide Wire lumen	10	Security key
3	Outer Sheat		10a Lock position
4	Stent		10b Extraction position
5	Stent stopper		10c Key after extraction
6	Coil	11	Pull back
7	Catheter body	12	Luer hub
8	Handle	13	Radiopaque Marker of the Sheat
		14	Radiopaque Markers of the Stent

English	2
Français	4
Deutsch	7
Italiano	9
Español	12
Dansk	14
Português	17
Suomi	19
Svenska	22
Nederlands	24
Ελληνικά	27
Türkçe	29
Česky	32
Eesti	34
Magyar	37
Latviski	39
Lietuvių kalba	42
Slovenčina	44
Slovenščina	47
Norsk	49
Македонски	52
Srpski	54
Hrvatski	57
Română	59
Polski	62
Български	64
Русский	67
简体中文	70
Bahasa Indonesia	72
Українська	75

ENGLISH

EASY FLYPE CARBOSTENT™ and EASY HIFLYPE CARBOSTENT™ iCarbofilm™ coated self-expanding peripheral stent systems

1. DEVICE DESCRIPTION

The Easy Flype Carbostent™ and Easy HiFlype Carbostent™ device family is composed of a self-expanding stent and a suitable system for correct placement and deployment of the former.

The **stent** is a flexible, self-expanding implantable device. It is made of Nickel and Titanium alloy (Nitinol) and coated with iCarbofilm™, a thin layer of carbon with a high density turbostratic structure substantially identical to that of the pyrolytic carbon used for mechanical cardiac valves discs.

Coating the substrate provides it with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Six radiopaque tantalum markers are applied to the stent (three at each end) to enable correct placement at the site of the target lesion.

The **delivery system** utilized is an over-the-wire type, ensuring that the stent can be positioned and released safely at the lesion. It is shown schematically in Figure 1 and consists of:

- a catheter body (7) 6F in diameter composed of:
 - a *retractable outer sheath*, the distal, semi-transparent section (3) of which contains the stent (4); the sheath has a radiopaque collar (13) near the tip of the catheter
 - an *inner lumen* comprising a coil (6) and an inner tube (2) within which a 0,035-inch guidewire can pass; the inner lumen extends along the catheter from a luer fitting (12) secured proximally to the handle and terminates distally at the tip (1) in radiopaque material, so as to facilitate the location and placement of the stent when viewed on an imaging system. A stent stopper (5) integral with the coil prevents axial displacement of the stent during retraction of the outer sheath.
- a handle (8) with an integrated delivery system, which allows the positioning of the catheter and release of the stent in situ. The handle has an ergonomic design that allows it to be grasped and used by the operator, even with one hand. It has a luer port (12) at the proximal end, which is used to flush the guidewire lumen, a safety key (10), a wheel (9) and a pull back tail (11).

The Manufacturer produces the peripheral stent and the deployment catheter system directly and implements all quality control procedures both during manufacture and on the finished products in accordance to the appropriate good manufacturing standards.

How supplied

The device is supplied sterile (using an ethylene oxide and CO₂ mixture), for single use, and individually packed in a pouch that must not be brought into a sterile area. Sterility is guaranteed while the package remains intact and until the expiry date printed on the packaging (EXPIRY DATE).

Contents

One (1) iCarbofilm™ coated self-expanding peripheral stent system.

Storage

The product should be stored as indicated on the packaging, in a cool and dry place, away from sunlight, at maximum temperature permitted up to 60 °C.

Do not use if the temperature indicator is black, as the free expansion diameter of the stent could be adversely affected.

2. INTENDED USE

The self-expanding peripheral stent system is intended for the treatment of the peripheral artery disease affecting the iliac artery, superficial femoral artery and the proximal third of popliteal artery.

3. INDICATIONS

The self-expanding peripheral stent system is indicated for the treatment of stenosed and occluded peripheral arteries (iliac artery, femoral artery and proximal third of the popliteal artery) in order to maintain vessel patency.

The self-expanding peripheral stent is also indicated in cases of suboptimal outcomes or dissections of the vessel wall resulting from angioplasty procedures.

4. CONTRAINDICATIONS

Use of the stent is contraindicated in the following situations:

- Lesions considered untreatable with PTA or other interventional techniques.
- Technically unsuccessful PTA procedure, for example due to the impossibility of accessing the stricture site with a delivery system.
- Patients with stenosis adjacent to an aneurysmal lesion of diameter at least twice the lumen of the native vessel.
- Lesions localized in the two distal thirds of the popliteal artery (or at the knee joint, generally considered).
- Clinical conditions, disorders or allergies that limit the use of anti-platelet and/ or anticoagulant therapy.
- Severe allergy to the contrast medium or drugs used during the procedure.

5. WARNINGS

- Evaluate the characteristics of the lesion to be treated and the specific physiopathology of the patient, with great care before making procedural choices.
- Considering the complexity of the procedure and the specific physiopathology of the patient, the physician should refer to updated literature for information regarding the risks and advantages of the various procedures before choosing the procedure to adopt.
- **In patients with known hypersensitivity or allergies to metal components**

of the stent, the device should be used only if potential benefits of the stent outweigh potential risks.

- There are no adequate or well-controlled studies in pregnant women for this product. The Easy Flype Carbostent™ and Easy HiFlype Carbostent™ devices should be used during pregnancy only if potential benefits of the stent outweigh potential risks.
- The device is intended for single use only. Do not reuse, reprocess or resterilize. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The device is designed to perform as a system. Do not use its components separately.
- Do not use the device if its packaging has been opened or damaged, incorrectly stored, or if the stated "Expiry date" has passed. In such cases product performances and sterility are not guaranteed.
- The device should be used only by physicians trained specifically to perform percutaneous transluminal angioplasty (PTA) and peripheral stent implantation.
- Follow the manufacturer's instructions when using accessories (introducer sheath, guidewire, haemostatic valve).
- The device should be guided under fluoroscopy and monitored with the use of radiographic equipment that produces high-quality images.
- If resistance is encountered at any stage of the procedure, do not force the system: withdraw the stent system and the entire delivery catheter, the guidewire and the introducer as a single unit. Applying excessive force and/ or incorrect handling may result in the stent being accidentally deployed or in damage to the placement catheter.
- Do not attempt to reposition a partially expanded stent. Attempted repositioning may cause severe vessel damage.
- At any time during operation, incorrect handling could cause dissection, perforation or rupture of the vessel.
- Do not advance or withdraw the stent system during deployment: it may cause serious damage to vessel walls.
- The stent cannot be recaptured or repositioned.
- The implantation of a self-expanding peripheral stent can cause dissection of the vessel distally and/or proximally to the stent and could also cause acute blockage of the vessel, making an additional operation necessary (surgical vascular replacement operation, further dilatation, deployment of additional stents or other procedures).
- Once the stent has been partly deployed, it cannot be recovered using the delivery system.
- Do not attempt to clean or re-sterilize devices that have been in contact with blood and organic tissues. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- The implantation of a stent through a bifurcation principal can obstruct or impede subsequent diagnostic or treatment procedures that may be envisaged.
- Extreme caution must be exercised in order to reduce the risk of plaque or other thrombotic material embolizing
- Where vessels are excessively tortuous and there is also extensive atheromatous plaque, the catheter may be difficult to advance. In such situations, incorrect handling could cause dissection, perforation or rupture of the vessel.
- Where two or more stents are needed, the materials of the single products must be of compatible composition.

6. PRECAUTIONS

- **Never use** the quick-release system at the start of stent expansion. Use it **only if necessary**, after releasing at least **half the length** of the stent.
- The device should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could produce mechanical alterations to it.
- The device's placement system is not intended for use with contrast medium injectors.
- Where stents are used to repair multiple lesions, implant them at the more distal sites first and work back toward the more proximally located, so as to avoid delivering catheters through stents already in place and risk displacing them.
- Where stents are implanted in tandem, the relative radiopaque markers should overlap. Overlap should not involve more than two stents.
- To avoid damaging the stent, use extreme care whenever passing the guidewire or a balloon catheter through stent that has just been expanded.

7. MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the manufacturer's Nitinol stent product lines are MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5T) or 3-Tesla (3T)
- Maximum spatial field gradient of 6,140 G/cm (61.40 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of:
 - 1.0 W/kg (Normal Operating Mode) at 1.5T
 - 1.0 W/kg (Normal Operating Mode) at 3T

RF Heating

Under the scan conditions defined above, the manufacturer's Nitinol stent product lines are expected to produce a maximum temperature rise of less than 6.8 °C after 15 minutes of continuous scanning at an MR system reported, whole body averaged specific absorption rate (SAR) of 1.0 W/kg.

	1.5 T	3 T
<i>Single Stent Configuration</i>		
<i>Highest temperature change</i>	≤6.3 °C	≤6.2 °C
<i>Overlapping Stent Configuration</i>		
<i>Highest temperature change</i>	≤6.8 °C	≤6.3 °C

MR Artifact

In non-clinical testing, the image artifact caused by the device extends approximately 0.7 cm from the manufacturer's Nitinol stents when imaged with gradient-echo and spin echo pulse sequences in a 3T MRI system.

8. DRUG REGIME

Studies in the clinical literature indicate the need to use appropriate anticoagulant and / or antiplatelet therapy with this stent system pre-operation, during the procedure and post-operation.

9. POTENTIAL ADVERSE EVENTS

The possible adverse events that may be associated with the implantation of a peripheral stent include risks related to percutaneous transluminal angioplasty and additional risks associated with the use of stents, as indicated below:

- Access-site infection, pain and haematoma
- Allergic reaction (to contrast medium, drugs used during the procedure or the materials composing the stent)
- Arterial spasm
- Arterial thrombosis
- Arteriovenous fistula
- Artery aneurysm or pseudoaneurysm
- Atherosclerotic distal embolization
- Bleeding
- Death
- Dissection, perforation or rupture of the stented vessel
- Drug (anti-coagulant/anti-platelet) and contrast medium-induced side effects/allergies
- Embolization (air, plaque, thrombotic material, stent or other)
- Hypo/hypertension
- Infection
- Kidney failure
- Lower limb ischemia (requiring surgical intervention)
- Restenosis of the artery in which the stent was implanted
- Stent occlusion
- Thrombosis
- Vessel occlusion or unresolved stenosis

10. DIRECTIONS FOR USE

Auxiliary material (not included in the box)

- Introducer with a **6F** diameter (**2.0 mm**) or greater
- Guidewires 0.035 inches in diameter (0.89 mm) or less
- Guiding catheter(s) with a **8F** diameter (inner diameter **2.7 mm**) or greater
- Syringes

MODELS

Each Easy Flype Carbostent™ and Easy HiFlype Carbostent™ device is identified by a model code and a lot number; available product codes are listed in *Table 1* and *Table 2*. The code is composed of the letters IC, followed by two letters identifying the type of device, two or three digits indicating the diameter of the expanded stent and a further two or three digits indicating the length of the stent, followed by the letter S or L depending on whether the catheter is 85 cm or 135 cm long.

The batch number allows traceability of all information regarding device manufacture and system control in the Quality Assurance archives of the Manufacturer.

To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included on each box; these labels can be attached to the implant patient's medical chart.

Selection of the device

Measure the length of the target lesion to establish the length of stent required. In order to appropriately cover the lesion, the deployed stent length should be at least 5 mm longer than the lesion itself on either side.

N.B.: percentage foreshortening of the stent is practically zero.

Measure the reference diameter of the vessel proximally and distally of the lesion to establish the appropriate model of stent.

To ensure a secure placement of the stent, the model selected must have a diameter at least one millimetre larger than the greater of the two reference diameters.

When selecting the model, refer to the information listed in *Table 1* and *Table 2*.

Preparing the device

- Open the box and take out the bag containing the stent and the delivery system. **Do not use if the temperature indicator is black, since the nominal diameter of the stent could be adversely affected.**
- After having examined the package for any damage, in a sterile area remove the support, on which the device is firmly held, from the bag. **Keep the support facing down.**
- Taking care to place the support on a flat surface in a sterile area, unfasten the band around the handle and remove the catheter from the protective coil. Check that the delivery catheter does not present kinks, bends or other damage. **Damage to the device may impair its performance.**
- Fill a 3 mL syringe with a physiological solution containing heparin, connect to the luer fitting for flushing (figure 1-12) and apply positive pressure. Continue flushing until the physiological solution emerges from the distal end of the catheter.
- Inspect the distal end of the catheter, making certain that the stent (figure 1-4) is covered completely by the outer sheath.

If the device is partially expanded or damaged in any way, do not use it.

Compatibility with accessories

Introducer sheaths of diameter 6F (2.0 mm) or greater can be used with the Easy Flype Carbostent™ and Easy HiFlype Carbostent™ devices. Guidewires of 0.035 inches (0.89 mm) or less can be used with the Easy Flype

Carbostent™ and Easy HiFlype Carbostent™ devices. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience. Should it be necessary to utilize a guiding catheter, systems are compatible with guiding catheters of diameter 8 F (2.7 mm) or greater. In this instance, make certain that the haemostatic valve is fully open before introducing the selected assembly.

Stent insertion

- Make sure the safety key (figure 1-10) of the handle is in the locked position (figure 1-10a).
- Wash the exposed portion of the guidewire with a physiological solution containing heparin, to remove traces of blood and contrast medium.
- Advance the system over the guidewire until the stent reaches the treatment site.

If unusual resistance is encountered at any time, do not force the system: withdraw the introducer, the guidewire and the device together as a single unit. Applying excessive force and/or incorrect handling may result in the stent being accidentally deployed or in damage to the placement catheter.


Advancing the device

- Advance the system to the point at which the stent is beyond the lesion.
- Draw the system back until the stent is in alignment with the lesion.
- Check that the part of the delivery system outside the body of the patient is freely extended and straight


Caution: if the delivery system is not properly extended, the stent could be deployed beyond the target lesion.

Stent implantation and expansion

- Firmly grasping the handle with one hand, check that the radiopaque markers on the stent are positioned proximally and distally in relation to the lesion.
- Make certain that the introducer sheath does not move during deployment of the stent.
- Turn the safety key (figure 1-10a) on the handle in the direction of the arrow (anticlockwise) to unlock it (figure 1-10b), then pull it out completely (figure 1-10c).
- Holding the handle in a fixed position with respect to the patient, check once again that the radiopaque markers are in the desired position with respect to the lesion; begin deploying the stent by turning the wheel (figure 1-9) with your thumb in the direction of the arrow, as shown in the figure:

- Under fluoroscopy, when the distal end of the stent is well-attached to the vessel wall, continue deployment by turning the wheel until the end of the stent is opened and it is completely implanted in the vessel walls.
- If rapid deployment is required, proceed as follows: make sure that there are no obstructions (e.g. female luer fitting screwed) in the distal section of the rear end of the handle, grasp the rear end of the device (figure 1-11) with your other hand and pull it, holding the front part of the handle with respect to the patient perfectly still, as shown in the figure:

- Check fluoroscopically that the stent has fully expanded.
- The stent can be expanded further, if necessary, using a PTA balloon catheter. In selecting the diameter of the balloon catheter and in performing the post-expansion, please refer to the diameter of the vessel that has been calculated by means of angiography. In any event, take particular care **not to expand the stent beyond its nominal diameter, stated on the label.**
- The final diameter of the stent must be suitable for the reference vessel diameter.

11. LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and packaged with the greatest care, using the most appropriate procedures the current state of technology allows. The safety standards integrated into the design and manufacture of the product guarantee its safe use under the aforementioned conditions and for its intended purposes in observation of the precautions listed above. These safety standards are meant to reduce as far as possible, but not eliminate completely, the risks related to the use of the product.

The product must be used only under a specialist physician's care, while taking into account any risks or side effects and complications that might arise from its

intended use, as mentioned in other sections of this instructions booklet. Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, the Manufacturer cannot be held responsible, either explicitly or implicitly, for the quality of end results subsequent to the use of the device or its effectiveness in resolving a patient's state of infirmity. The end results, in terms of both the clinical status of the patient and the functionality and the lifetime of the device, depend on many factors outside the manufacturer's control, among which are the patient's conditions, the surgical procedure of implantation and application, and the handling of the device after being removed from its package.

In light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. Under such circumstances, the customer shall deliver the device to The Manufacturer, which reserves the right to examine the allegedly faulty device and establish whether it actually presents manufacturing defects. The warranty consists exclusively in the replacement of the defective device with another device of the same type or equivalent, from the same Manufacturer.

The warranty applies under the condition that the device is returned correctly packed to the Manufacturer, and that it is accompanied by a written, detailed report describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient.

On replacement of the device, the Manufacturer shall reimburse the buyer the expenses incurred for return of the defective device.

The Manufacturer declines any responsibility for cases of negligence to observe the methods of use and the precautions stated in this instructions booklet and for cases of use of the device after the use before date printed on the package. Furthermore, the Manufacturer declines any responsibility related to the consequences arising from the treatment choices and the methods of use or the application of the device; the Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, subsequent to the application of the device, or the choice of implantation technique used by the operator.

The agents and the representatives of the Manufacturer are not authorized to modify any of the conditions of this warranty or to assume any further obligations or to offer any guarantees related to this product beyond the terms stated above.

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