Invader™ PTCA

PTCA BALLOON DILATATION CATHETER

EN	ENGLISH.....	5	KO	한국어.....	41
TR	TÜRKÇE	6	SL	SLOVENŠČINA.....	42
FR	FRANÇAIS.....	7	FI	SUOMI.....	44
DE	DEUTSCH	9	ET	EESTI.....	45
NL	NEDERLANDS.....	10	LV	LATVIEŠU	46
ES	ESPAÑOL.....	11	LT	LIETUVIŲ.....	48
IT	ITALIANO	13	NO	NORSK.....	49
PT	PORTUGUÊS	14			
HU	MAGYAR.....	15			
HR	HRVATSKI	17			
BS	BOSANSKI.....	18			
CS	ČEŠTINA	19			
DA	DANSK.....	21			
SR	SRPSKI.....	22			
SV	SVENSKA.....	23			
RU	РУССКИЙ.....	25			
PL	POLSKI.....	26			
RO	ROMÂNĂ	28			
MK	МАКЕДОНСКИ.....	29			
BG	БЪЛГАРСКИ	30			
EL	ΕΛΛΗΝΙΚΑ.....	32			
SK	SLOVENČINA	33			
ZH	简体中文	34			
JA	日本語.....	35			
ID	BAHASA INDONESIA	36			
AR	العربية	38			
FA	فارسی	39			
UK	УКРАЇНСЬКА.....	40			

INSTRUCTIONS FOR USE

DESCRIPTION

STORAGE

INTENDED USE

INDICATIONS

CONTRAINDICATIONS

WARNINGS / PRECAUTIONS

POTENTIAL ADVERSE EFFECTS

INSTRUCTIONS FOR USE

Auxiliary material (not included in the package)

- Selecting the right size
- Preparation
- Insertion
- Removal

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ENGLISH

INVADER™ PTCA BALLOON DILATATION CATHETER HYDROPHILIC-COATED PTCA BALLOON DILATATION CATHETER

DESCRIPTION

INVADER™ PTCA Balloon Dilatation Catheter is a single use, Ethylene Oxide (EO) sterilized, monorail rapid exchange (RX) coronary angioplasty balloon dilatation catheter. The ALVIMEDICA INVADER™ PTCA has been designed to dilate the stenotic atherosclerotic lesions in coronary arteries or bypass grafts.

The dilation part of the catheter is the balloon near the distal tip. The catheter is hydrophilic coated on the distal shaft excluding the balloon. Radiopaque marker bands are located on both the proximal and distal shoulders of the balloon. Two depth markers, on the proximal shaft aid in determining when the balloon is exiting the guiding catheter, in the cases of brachial or femoral approach, respectively. A separate lumen on the catheter shaft is intended for use as a guidewire lumen beginning at approximately 27 cm from the distal tip. The PTFE coated proximal end of the catheter is used as the balloon inflating port. The balloon is inflated by injecting a contrast material from this end. The balloon material can be inflated to a given size at a given pressure. The size matrix and device features are as below:

Table 1. The size matrix and device features of INVADER™ PTCA Balloon Dilatation Catheter

Balloon Outer Diameter (mm)	Balloon Length (mm)	Minimum Guiding Catheter Compatibility ID	Nominal Pressure (bar)	Rated Burst Pressure (bar)
1.25	6, 10, 12, 15, 20	5F (1.4 mm / 0.056")	8	16
1.50	6, 10, 12, 15, 20	5F (1.4 mm / 0.056")	8	16
2.00	6, 10, 12, 15, 20, 25	5F (1.4 mm / 0.056")	8	16
2.25	6, 10, 12, 15, 20, 25	5F (1.4 mm / 0.056")	8	16
2.50	6, 10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
2.75	6, 10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
3.00	6, 10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
3.25	10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
3.50	10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
3.75	10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
4.00	10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
4.50	10, 12, 15, 20, 25, 30, 34	5F (1.4 mm / 0.056")	8	16
5.00	10, 12, 15, 20	5F (1.4 mm / 0.056")	8	16

STORAGE

Store in a cool, dry place, away from sunlight, at a temperature between 10°C and 40°C.

INTENDED USE

The INVADER™ PTCA Balloon Dilatation Catheter is intended to be used for dilating stenoses in the coronary artery or bypass graft so as to increase myocardial perfusion.

- Patients should be eligible for coronary bypass surgery.
- It is indicated for patients with single-artery non-calcified atherosclerotic lesions that can be dilated using a PTCA catheter.
- This procedure can also be indicated in certain patients who have multi-artery disease, and in patients who have undergone aorta-coronary bypass surgery but still have:
 - > recurrent symptoms;
 - > progressive coronary artery disease;
 - > stenosis or obstruction in bypass grafts.

INDICATIONS

Patients who are eligible for percutaneous coronary intervention (PCI).

Product should not be used in the pregnant or nursing women and pediatric patients.

CONTRAINDICATIONS

The INVADER™ PTCA Balloon Dilatation Catheter is contraindicated in patients:

- who are not eligible for coronary bypass surgery;
- who have fully obstructed coronary arteries;
- who have diffuse lesions;
- who have severe stenosis of the left main coronary artery.

WARNINGS / PRECAUTIONS

- The INVADER™ PTCA Balloon Dilatation Catheter has been designed for single use only, reuse is not recommended. Do not resterilize it.
- Keep the catheter in a cool, dry and dark place.
- Do not use the catheter after the expiry date printed on the packaging.
- Use diluted contrast material only.
- Do not use air or any other gas to inflate the balloon.
- Check the packaging for any damage.
- Do not exceed the rated burst pressure as indicated in the Instructions For Use when inflating the balloon.
- Use of an inflation device with an incorporated pressure gauge is recommended.
- The diameter of the inflated balloon should not be exceeded at the points just proximal and distal to the stenosis.
- The INVADER™ PTCA Balloon Dilatation Catheter should only be used by experienced physicians who have been trained in PTCA operations.
- Give appropriate anticoagulation and vasodilatation therapy prior to catheterization.
- The PTCA operation should only be performed in medical centers capable of carrying out emergency coronary bypass surgery in case of severe complications.
- Do not tighten the hemostatic adapter in the Y-connector, as this can compress the shaft, thus impeding the inflation and deflation of the balloon.
- All procedures performed once the catheter has been introduced into the body should be carried out under quality fluoroscopy. Never pull or push the catheter unless the balloon has been fully deflated under vacuum. If any resistance is encountered during the procedure, simply stop and try to identify the cause and then gradually advance the balloon catheter to prevent the occurrence of kinking on the proximal shaft. If you fail to identify the cause, remove the entire system.
- Do not use the contrast materials Ethiodol or Lipiodol.
- Do not expose the insertion system to organic solvents (i.e. alcohol, etc.)
- If a serious incident related to device occurs, immediately report the incident to Alvimedica and applicable competent authority or regulatory body.

The product might be a potential biohazard, after use. The product must be handled and disposed of in accordance with accepted medical practice, applicable local laws and regulations. Alvimedica is not responsible for the handling and disposal of the product after use.

Summary of safety and clinical performance

The Summary of Safety and Clinical Performance (SSCP) for this device can be found at alvimedica.com and in EUDAMED database when available.

POTENTIAL ADVERSE EFFECTS

- Dissection of the coronary artery
- Tearing, perforation of, or damage to the coronary artery
- Complete obstruction of the coronary artery or bypass graft
- Thrombosis of the coronary artery
- Unstable angina
- Acute myocardial infarction
- Restenosis of the dilated artery
- Spasm of the coronary artery
- Arrhythmias, including ventricular fibrillation
- Hemorrhage and hematoma
- Drug reactions, allergic reaction to contrast material
- Hypertension - hypotension
- Infection

- Arteriovenous fistula
- Embolism
- Death
- Urgent coronary artery bypass graft surgery
- Renal failure

INSTRUCTIONS FOR USE

Auxiliary material (not included in the package)

- Guiding catheter(s) with a 5F diameter (1.422 mm) or larger
- Syringes
- Guidewires 0.014 inches (0.356 mm) in diameter or less
- Inflation device
- Introducer sheath

a. Selecting the right size

The diameter of the balloon should be selected so that, when it is inflated at nominal pressure, it will not exceed the coronary artery diameter at the points just proximal and distal to the stenosis.

b. Preparation

- Prepare an inflation solution of a 1:1 mixture of contrast material and sterile saline solution.
- Fill a 20 cc syringe with approximately 10 cc of the saline solution.
- Place a needle on the syringe and then carefully insert the needle through the distal tip of the catheter and flush the guidewire lumen.
- Attach a stopcock to the hub of the INVADER™ PTCA Balloon Dilatation Catheter.
- Attach the syringe to the stopcock and then hold the syringe vertically so that its plunger is at the top. Pull the plunger and purge the balloon of air.
- Apply a vacuum using the syringe for 15-20 seconds and make sure that no air bubbles pass through the diluted contrast material.
- Carefully stop vacuuming.
- Repeat the process if required.
- Exert negative pressure and then switch off the stopcock.
- Make sure that no air enters the system; attach the inflation device to the catheter. Switch on the stopcock, inflate the balloon to nominal pressure and then check the integrity of the INVADER™ PTCA Balloon Dilatation Catheter.
- Exert negative pressure and then switch off the stopcock.

c. Insertion

- Using a hemostatic valve, attach the Y-connector to the guiding catheter already inserted in the femoral artery.
- Pass the guidewire (maximum 0.014" (0.356 mm)) through the guiding catheter, advance it and then place it as required.
- Gently tighten the hemostatic valve of the Y-connector on to the guidewire so as to prevent unintended movement of the guidewire.
- Insert the rear tip of the guidewire through the distal tip of the catheter and attach it to the catheter. The rear tip of the guidewire will come out at approximately 27 cm proximal to the catheter's distal tip.

WARNING: Kissing balloon: bench tests has shown that two INVADER™ PTCA Balloon Dilatation Catheter with a diameter ≤ 4mm can be introduced in 6F (or larger) guide catheter.

- Loosen the hemostatic valve of the Y-connector and advance the INVADER™ PTCA Balloon Dilatation Catheter towards the distal part of the guiding catheter.

WARNING: always advance the INVADER™ PTCA Balloon Dilatation Catheter when it is deflated and is on the guidewire.

- Do not tighten the hemostatic valve of the Y-connector too much, as this may prevent the passage of the contrast material through the balloon inflation lumen, thus resulting in the prolongation of the deflation/inflation time of the balloon.
- Using standard angioplasty techniques, continue the procedure under fluoroscopy. Radiopaque markers can help to better position the balloon in the stenosis.

Note: It is recommended that the guidewire and/or balloon catheter remain in the lesion until the dilatation is complete. Variations in the viscosity of the contrast material may affect the deflation/inflation time.

TÜRKÇE

d. Removal

- Loosen the hemostatic valve of the Y-connector.
- Hold the hemostatic valve and the guidewire with one hand and the shaft of the catheter with the other.
- Whilst keeping the guidewire immobile in order to maintain its position in the coronary artery, remove the dilatation catheter from the guiding catheter.

Note: During removal, check the position of the guidewire by means of fluoroscopy.

- Switch on the hemostatic valve of the Y-connector.
- If necessary, prepare another INVADER™ PTCA Balloon Dilatation Catheter and repeat the dilatation process.

RESPONSIBILITY AND WARRANTY

Alvimedica guarantees that each and every component of this product has been manufactured, packed, tested and sterilized without any defect to its workmanship or material. Each product has been tested prior to packaging. Alvimedica shall exchange any product with manufacturing or packaging defect(s) with a new one upon the decision of a competent binding judicial authority. Due to biological variations among individuals, no product is 100 % effective in every case. Therefore, and since Alvimedica has no control over conditions under which the product is used after sales, the selection of patients and the methods of application, Alvimedica disclaims all warranties with respect to the product, expressed or implied including, without limitation, those of merchantability, and fitness for a particular purpose. Alvimedica is not directly or indirectly responsible for any injury or damage to or loss of any person resulting from the use of the product, nor is Alvimedica responsible or liable, directly or indirectly, for any injury, damage or loss that may result from re-use or re-sterilization.

LIMITATION OF LIABILITY

In no event shall Alvimedica be liable for any special, consequential, indirect, collateral, incidental or punitive damages or lost profits or failure to realize expected savings or other commercial or economic loss of any kind, arising out of any breach of this IFU including breach of any precautions, warnings and/or any recommendation mentioned herein.

INTELLECTUAL PROPERTY PROTECTION

All trademarks, logos, inventions, know-how, technology, proprietary information and other intellectual property rights of this product belong solely to Alvimedica. Alvimedica reserves all rights to demand, sue, seek relief before judicial authorities in case of an infringement regarding any unauthorized use or misuse of intellectual property or reverse engineering related to the product.

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