Alvision™

INTERVENTIONAL CARDIOLOGY DIAGNOSTIC CATHETER

EN	ENGLISH	2
TR	TÜRKÇE.....	3
FR	FRANÇAIS	4
DE	DEUTSCH	6
NL	NEDERLANDS.....	7
ES	ESPAÑOL.....	8
IT	ITALIANO	10
PT	PORTUGUÊS.....	11
HR	HRVATSKI	13
HU	MAGYAR.....	14
BS	BOSANSKI.....	16
CS	ČEŠTINA	17
DA	DANISH	19
SR	SRPSKI.....	20
SV	SWENSKA	22
RU	PYСCKИЙ.....	23
PL	POLSKI	24
RO	ROMÂNĂ	26
MK	МАКЕДОНСКИ.....	27
BG	БЪЛГАРСКИ	29
EL	Ελληνικά	30
SL	SLOVENŠČINA.....	32
SK	SLOVENČINA.....	33
ZH	中文	35
JA	日本語	36
AR	العربية	39
FA	فارسی	41

Rx only	Prescription only-device restricted to use by or on order of a physician.
	Caution.
	Consult instructions for use.
STERILE EO	Sterilized with Ethylene Oxide.
	Keep between 10-40°C.
	Keep away from sunlight.
	Keep dry.
	Do not reuse.
	Do not resterilize.
	Do not use if package is damaged.
	Non-Pyrogenic
	Contents
	Inner Diameter
	Length
PSI / kPa	Max. Pressure
	Max. Guidewire
	Date of Manufacture
	Use by
REF	Catalogue Number
LOT	Batch Code
	Manufacturer
EC REP	Authorised Representative in the European Community

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ENGLISH

Instructions for Use – Alvision™ Interventional Cardiology Diagnostic Catheters

Description

Alvision™ Interventional Cardiology Diagnostic Catheters incorporate an atraumatic tip. These catheters are available in various sizes and configurations. They have been sterilized with ethylene oxide gas, non-pyrogenic. They are intended for single use only.

Specifications

The specifications can be found on the label of each catheter pouch.

Indications

Alvision™ Interventional Cardiology Diagnostic Catheters are intended for use in the delivery of radio-opaque media to selected sites in the coronary vascular system.

Contraindications

Passing catheters through synthetic vascular grafts should be avoided.

Warnings / Precautions

- Diagnostic catheters are invasive devices that pose a potential risk of injury or death. Accordingly, angiographic catheters should only be used by trained, experienced physicians who are aware of the potential complications and risks of angiography.
- Extreme care must be taken to avoid causing damage to the vessel walls through which the catheter passes.
- This catheter may occlude smaller vessels.
- Take care to avoid completely blocking the blood flow.
- Do not reuse or reshape the catheter.
- The maximum permitted diameter of the guidewire is specified on the pouch. Guidewire use has been associated with thrombus formation.
- Alvision™ Interventional Cardiology Diagnostic Catheters are intended for single use only and should be disposed of after one procedure.
- Do not allow catheters to come into contact with organic solvents.
- Do not pressurize catheters to in excess of the rated maximum pressure marked on the pouch label.
- A thorough inspection should be made to verify whether the size, shape and condition of the catheter are suitable for the procedure concerned.
- Angiography should only be performed by an experienced physician. Catheters are supplied sterile and non-pyrogenic.
- Do not use if the package has been opened or damaged.
- As with any surgical procedure, it is recommended to apply "Universal Precautions".
- Avoid damaging the tip of the catheter during removal from the packaging.
- Removing too rapidly may result in damage to the catheter tip.
- Take care to prevent or reduce clotting when using a catheter.
- Whilst in the vascular system, the catheter should always be filled with either flushing solution or contrast medium.
- The use of systemic heparinization should be considered.
- These catheters are designed for use in the injection of contrast media as well as to facilitate physiological measurement in conjunction with routine diagnostic procedures.
- Do not use after the "Use By" date specified on the label.
- Do not expose catheters to strong solvents or ionizing radiation.
- Do not store in excessive heat.
- Avoid prolonged exposure to UV light.
- Do not re-use, re-process or re-sterilize. This will affect the product performance. The product may not be able to withstand the stresses of re-use. Cleaning or re-processing of the catheter components exposed to biological materials may result in adverse device effect when used.

The product might be a potential biohazard, after use. Handling and dispose of products in accordance with accepted medical practice, applicable local laws and regulations must be done and is not in the responsibility of Alvimedica.

Potential Adverse Events

The following adverse reactions have been reported to be associated with the use of diagnostic catheters: thrombus formation/emboli, arterial wall damage, plaque dislodgement, myocardial infarction, cardiac arrhythmia, stroke and death.

Physicians should be familiar with the current literature concerning the complications associated with angiography.

Possible adverse events include:

- Air embolism
- Hematoma
- Infection
- Perforation of the heart or vessel wall

Instruction for Use

- Gently remove the catheter hub first from the pouch so as to prevent damaging the tip.
- The catheter should be introduced into the blood vessel using the physician's preferred vascular entry technique.
- Guidewires should be used to straighten the catheters and prevent kinking during introduction.
- Thin-walled catheters should be handled with extreme care in order to prevent kinking.

TERMS OF WARRANTY

Alvimedica guarantees that each and every component of this product has been manufactured, packed, tested and sterilized without any defect to its workmanship or material. Each product has been tested prior to packaging. Alvimedica shall exchange any product with manufacturing or packaging defect(s) with a new one upon a decision of a competent binding judicial authority. Due to biological variations among individuals, no product is 100 % effective in every case. Therefore, and since Alvimedica has no control over conditions under which the product used after sales, the selection of patients and the methods of application; Alvimedica disclaims all warranties with respect to the product, expressed or implied including, without limitation, those of merchantability, fitness for a particular purpose. Alvimedica is not directly or indirectly responsible for any injury or damage to or loss of any person resulting from the use of the product, nor is Alvimedica responsible or liable, directly or indirectly, for any injury, damage or loss that may result from re-use or re-sterilization.

LIMITATIONS OF LIABILITY

In no event shall Alvimedica be liable for any special, consequential, indirect, collateral, incidental or punitive damages or lost profits or failure to realize expected savings or other commercial or economic loss of any kind, arising out of any breach of this IFU including breach of any precautions, warnings and/or any recommendation mentioned herein.

INTELLECTUAL PROPERTY PROTECTION

All trademarks, logos, inventions, know-how, technology, proprietary information and other intellectual property rights of this product solely belongs to Alvimedica. Alvimedica reserves all the rights to demand, sue, seek relief before judicial authorities in case of an infringement regarding any unauthorized use or misuse of intellectual property or reverse engineering related to the product.

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