# Alviguide™ Blue+

INTERVENTIONAL CARDIOLOGY GUIDING CATHETER

EN	ENGLISH .....	2
TR	TÜRKÇE.....	3
FR	FRANÇAIS .....	5
DE	DEUTSCH .....	7
NL	NEDERLANDS.....	9
ES	ESPAÑOL.....	11
IT	ITALIANO .....	13
PT	PORTUGUÊS.....	14
HR	HRVATSKI .....	16
HU	MAGYAR.....	18
BS	BOSANSKI .....	20
CS	ČEŠTINA.....	22
NO	NORSK .....	23
DA	DANSK.....	25
SR	SRPSKI .....	27
SV	SVENSKA.....	29
RU	РУССКИЙ .....	31
PL	POLSKI .....	33
RO	ROMÂNĂ .....	34
MK	МАКЕДОНСКИ .....	36
BG	БЪЛГАРСКИ .....	38
EL	ΕΛΛΗΝΙΚΑ.....	40
SL	SLOVENŠČINA.....	42
SK	SLOVENČINA.....	44
ET	EESTI.....	45
FI	SUOMI .....	47
LT	LIETUVIŲ .....	49
LV	LATVISKI.....	50
AR	العربية .....	53
FA	فارسی .....	54
ZH	中文 .....	56
JA	日本語 .....	57

**Rx only**

Prescription only-device restricted to use by or on order of a physician.

Caution.

Consult instructions for use.

**STERILE EO**

Sterilized with Ethylene Oxide.

Keep between 10-40°C.

Keep away from sunlight.

Keep dry.

Do not reuse.

Do not resterilize.

Do not use if package is damaged.

Non-Pyrogenic

Contents

Inner Diameter

Length

Date of Manufacture

Use by

**REF**

Catalogue Number

**LOT**

Batch Number

Manufacturer

**EC REP**

Authorised Representative in the European Community

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### ALVIGUIDE™ BLUE PLUS INTERVENTIONAL CARDIOLOGY GUIDING CATHETER

#### DEVICE DESCRIPTION

The Alviguides™ Blue Plus Interventional Cardiology Guiding Catheter is a single-use product that provides optimal back-up support and control for various anatomies with its large lumen. The catheter is designed to provide a pathway through which therapeutic and diagnostic devices are introduced. It is available in various sizes (5F, 6F, 7F, 8F) with the catheter length of 100 cm.

#### HOW SUPPLIED

The Alviguides™ Blue Plus Interventional Cardiology Guiding Catheter is supplied with a single Instructions For Use (IFU) in each packaging box.

**Sterile:** This device is sterilized with Ethylene Oxide.

**Single use:** This device is for single use only. Do not re-sterilize. Do not use if the package is open or damaged.

#### INDICATIONS

The Alviguides™ Blue Plus Interventional Cardiology Guiding Catheter is intended for use in the intravascular introduction of interventional/diagnostic devices into the coronary vascular system.

#### CONTRAINDICATIONS

Passing catheters through synthetic vascular grafts should be avoided.

#### WARNINGS/PRECAUTIONS

- The product should only be used by experienced physicians who have been adequately trained in percutaneous coronary intervention techniques.
- The procedure should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed in the event of potentially injurious or life-threatening complications.
- Do not use the product when integrity of the sterile packaging has been compromised or if the device appears damaged.
- Before using the product, read all instructions for use thoroughly prior to procedure.
- Do not use the product after the use-by date.
- Do not re-use, re-process or re-sterilize. This will affect the product performance. The product may not be able to withstand the stresses of re-use. Cleaning or re-processing of the catheter components exposed to biological materials may result in adverse device effect during reuse.
- Due to the size and relative stiffness of guiding catheters, extreme care must be taken to avoid causing damage to the vessel walls through which this catheter passes. Inspect the guiding catheter prior to use in order to verify that its size, shape and condition are suitable for the procedure concerned.
- Do not use the contrast materials Ethiodol or Lipiodol.
- Do not expose to organic solvents.
- It is recommended to remove the peel-off label for traceability purposes before disposal of the packaging.

If strong resistance is met during manipulation, discontinue the procedure and determine the cause of the resistance before proceeding. If the cause of the resistance cannot be determined, withdraw the catheter.

- Advancement, manipulation and withdrawal of the guiding catheter should always be performed under fluoroscopic visualization.
- Sideholes should not be placed in the shaft of the guiding catheter by the user. Puncturing the shaft of the guiding catheter with hospital instruments may lead to thrombogenesis or failure of shaft integrity.
- Prior to use, examine the catheter to verify proper function and ensure that its size and shape are suitable to the specific procedure for which it is to be used.
- Appropriate anticoagulation of the patient is indicated with the use of this device.
- The product might be a potential biohazard after use. Handling and disposal of products must take place in accordance with accepted medical practice, applicable local laws and regulations and is not in the responsibility of Alvimedica.

#### POTENTIAL ADVERSE EVENTS/COMPLICATIONS

- Additional intervention (percutaneous/surgical)
- Air embolism
- Arrhythmia
- Damage to healthy intima
- Dissection
- Embolism
- Hemorrhage/Hematoma

- Infection
- Myocardial infarction
- Myocardial ischemia
- Stroke
- Thrombosis
- Total occlusion of coronary artery
- Vessel spasm

#### INSTRUCTIONS FOR USE

1. Remove the guiding catheter from its packaging.
2. Upon removal from its packaging, inspect the guiding catheter to verify whether it has been damaged.

**Note:** Do not use a guiding catheter that has been damaged in any way. If damage is detected, replace with an undamaged guiding catheter.

3. Prior to use, flush the lumen of the guiding catheter with a heparinized saline solution.
4. Introduce the guiding catheter into the vasculature via the catheter introducer, and over a guidewire using the vessel entry technique of choice.
5. Advance the catheter and guidewire together to the selected vascular site under fluoroscopic visualization.
6. Set up a continuous heparinized saline flush attached to the guiding catheter hub. It is recommended that a continuous heparinized saline flush be maintained between the guiding catheter and any intraluminal device that passes through it.
7. To infuse the contrast media or to introduce other intravascular devices the diagnostic guidewire has to be removed.

#### TERMS OF WARRANTY

Alvimedica warrants that reasonable care has been taken in the design and manufacture of this instrument. Due to biological variations among individuals, no product is 100% effective in every case. Therefore, and since Alvimedica has no control over conditions under which the product is used after sales, the selection of patients and the methods of application, Alvimedica disclaims all warranties with respect to the product, expressed or implied including, without limitation, those of merchantability, fitness for a particular purpose.

Alvimedica is not directly or indirectly responsible for any injury or damage to or loss of any person resulting from the misuse of the product, nor is Alvimedica responsible or liable, directly or indirectly, for any injury, damage or loss that may result from re-use or re-sterilization.

#### INTELLECTUAL PROPERTY PROTECTION

All trademarks, logos, inventions, know-how, technology, proprietary information and other intellectual property rights of this product belong solely to Alvimedica. Alvimedica reserves all the rights to demand, sue or seek relief before judicial authorities in case of an infringement regarding any unauthorized use or misuse of intellectual property or reverse engineering related the product.

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