
CiD
Carbostent® Implantable Devices

INTERVENTIONAL CARDIOLOGY

CARBOFILM™ TECHNOLOGY FOR A BIO INDUCER SURFACE

AVANT GARDE

COBALT CHROME CHRONO CARBOSTENT

Carbofilm Coated Coronary Stent on RX balloon catheter
Stent coronaire recouvert de Carbofilm sur cathéter à ballonnet à échange rapide
Karbonbeschichteter Koronarstent auf einem RX Ballon-Katheter
Stent coronarico rivestito Carbofilm su catetere a palloncino a scambio rapido
Stent coronario revestido de Carbofilm en catéter de balón de intercambio rápido
Stent coronário revestido de Carbofilm em cateter de balão de troca rápida
Met carbonfilm beklede coronaire stent op balloncatheter met snelle uitwisseling
Carbofilm-coated koronarstent på RX ballonkateter
Koronarstent med yttsikt av Carbofilm på ballongkateter för snabbbyte
Carbofilm-päällysteinen pikavaihdeettavalle pallokatetrille asennettu koronaaristent
Στεφανιαίο Stent επικαλυμμένο με Carbofilm σε καθετήρα μπαλονιού ταχείας αλλαγής
Adet RX balon kateter üzerinde Karbofilm Kaplı Koroner Stent
Jeden koronární stent potažený karbofilmem na balonkovém katetru RX
Süsinikkilega kaetud koronaarstent RX balloonkatetritil
Karbofilm bevonatú koronária sztent RX ballon katéteren

„Carbofilm“ padengtas vainikinių arterijų stentas ant „RX“ balioninio kateterio
Koronární stents ar Carbofilm pārklājumu, uzmontēts uz RX balonkatetra
Koronarstent dekket av karbonfilm på RX ballongkateter
Uhlíkovým povlakom potiahnutý koronárny stent na RX balónikovom katétri
Koronarni stent s Carbofilm prevleko na hitro zamenljivem balonskem katetru
Коронарен стент прекриен со Carbofilm на RX балон катетер
Karbofilmom obloženi koronarni stent na RX balon kateteru
Karbofilmom obloženi koronarni stent na RX balon kateteru
Stent coronar acoperit cu carbofilm pe cateterul cu balon RX
Wieńcowy Stent Pokryty Carbofilmem na RX baloniku cewnika
Коронарен стент с покритие Carbofilm върху RX балонен катетър.
Коронарний стент з покриттям Carbofilm на балонному катетері швидкої заміни (RX)
Коронарний стент с покриттям Carbofilm на баллонном катетере быстрой замены (RX)
RX球囊导管Carbofilm涂层冠状动脉支架

IC0710.137.105 Issued 202001 Rev.9

CiD
Carbostent & Implantable Devices

Manufactured by: **CiD** S.p.A.
Strada per Crescentino, sn - 13040 Saluggia (VC), Italy
CiD S.p.A. member of **Alvimedica** group
www.alvimedica.com

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DESCRIPTION

The "AVANTGARDE CHRONO™ CARBOSTENT" device consists of a coronary stent firmly held on the distal end of a semicompliant balloon catheter

The **coronary stent** is a flexible implantable device that can be expanded using a PTCA catheter.

The stent is made of cobalt chromium alloy coated with Carbofilm™, a thin turbostratic carbon film with a high-density crystalline structure substantially identical to that of the pyrolytic carbon used for mechanical cardiac valve discs.

The coating provides the substrate with the bio- and haemocompatible characteristics of pyrolytic carbon, without affecting the physical and structural properties of the substrate itself.

Two radio-opaque platinum markers are located at either end of the stent and enable accurate positioning of the device over the lesion to be treated.

The rapid-exchange type **balloon catheter** provides a safe means for conveying the coronary stent to the lesion to be treated.

The distal portion of the catheter consists of two lumens: one is used for balloon inflation and deflation, the other for the passage of the guidewire.

Two radio-opaque markers, located outside the usable length of the balloon enable accurate placement across the stenosis.

The proximal portion of the catheter, with stainless steel hypotube, forms the lumen for balloon inflation and deflation.

Two depth indicators placed 90 and 100 cm from the distal end for femoral and brachial approach, respectively, aid in determining when the balloon is exiting the guiding catheter. The proximal end of the catheter has a female Luer lock for attachment to an inflation device. The Manufacturer directly manufactures the "AVANTGARDE CHRONO™ CARBOSTENT" device and carries out all quality controls both during production and on the finished product, according to the norms of Good Manufacturing Practice.

INTENDED USE

Use of the coronary stent is indicated in the treatment of coronary occlusions in order to maintain patency of the vessel.

INDICATIONS

The stent is indicated for use in the following cases:

- Elective treatment of obstructive lesions of native coronaries and aortocoronary bypasses in patients having the requisites for percutaneous transluminary coronary angioplasty (PTCA)
- When short- and long-term results of PTCA procedures are inadequate
- Primary angioplasty or intervention during acute myocardial infarction.

CONTRAINDICATIONS

The use of the stent is contraindicated in the following situations:

- Patients who do not have the requirements for an aortocoronary bypass operation
- Pregnant women
- Disorders that limit the use of antiaggregant and/or anticoagulant therapy
- Severe allergy to contrast media
- Stenosis of the unprotected vessel
- Ejection fractions $\leq 30\%$
- Lesions considered untreatable with PTCA or other operating techniques
- Resistant lesions which cannot be dilated
- Lesions on a vessel with a reference diameter < 2.25 mm

MODEL

Each "AVANTGARDE CHRONO™ CARBOSTENT" device is identified by a model code and a lot number; available models are listed in Table 1. The model code consists of the letters AG, which specify the name of the device, followed by two or three digits indicating the nominal expansion diameter of the stent, and by two other digits indicating the stent length. The lot number allows traceability of all information regarding the device manufacturing process and system control in the Quality Assurance archives of the Manufacturer. To facilitate traceability of the device at the user end, the product code is printed on the adhesive labels included in each box; these labels can be attached to the patient's medical chart.

PACKAGING

The "AVANTGARDE CHRONO™ CARBOSTENT" is supplied sterile and individually packed in a pouch that **must not be placed in a sterile field**.

The sterilisation made by the manufacturer uses an ethylene oxide and CO₂ mixture. Sterility is guaranteed while the package remains intact and until the expiry date printed on the packaging (EXPIRY DATE).

WARNINGS

- The device is intended for single use only. Do not reuse, reprocess or resterilise. This could cause contamination risks of the device and/or infections to the patient, inflammation and patient-to-patient transmission of infective diseases.
- The "AVANTGARDE CHRONO™ CARBOSTENT" device is designed to perform as a complete system. Do not use its components separately.
- The "AVANTGARDE CHRONO™ CARBOSTENT" device is also indicated for use in association with PTCA. The positioning catheter is not a coronary dilation catheter; it should be used for stent use only.
- Do not use the "AVANTGARDE CHRONO™ CARBOSTENT" device if its packaging has been opened or damaged or if the period of sterility has expired. In such cases product sterility is not guaranteed.
- Do not remove the stent from the release balloon, as this could damage the stent and/or cause its embolization.
- The "AVANTGARDE CHRONO™ CARBOSTENT" should be carefully handled so as to avoid any contact with metallic or abrasive instruments that could damage the highly polished surfaces or produce mechanical impairment.

- Do not use the "AVANTGARDE CHRONO™ CARBOSTENT" device if, during handling, due to rotations or forced insertion, the proximal portion of the catheter has kinks or bends; in such cases do not try to straighten the catheter.
- The "AVANTGARDE CHRONO™ CARBOSTENT" device should be used by physicians specifically trained in performing PTCA and in implanting coronary stents.
- A cardiac surgery team should be on stand-by.
- The system should be handled under fluoroscopic control using X-ray equipment that produces high-quality images.
- The characteristics of the lesion to be treated and the specific physiopathology of the patient should be assessed very carefully before taking procedural decisions. The physician in charge of supervising the system may decide to carry out preliminary dilation for release of the stent if there are lesions that permit this approach.
- The stent must be implanted in the target lesion using the positioning catheter on which it is mounted.
- Do not pressurise the system until the stent is positioned across the lesion to be treated.
- Do not retract the catheter at the end of the procedure until its balloon is fully deflated.
- Do not attempt to reposition a partially expanded stent. This action may cause severe vessel damage.
- Do not attempt to clean or resterilize devices that have been in contact with blood and organic tissue. Used devices are to be disposed of as dangerous medical waste with risk of infection.
- Follow the manufacturer's instructions for the use of accessories (guide catheter, guidewire, haemostatic valve).
- If unusual resistance is encountered at any time during the insertion procedure, do not force the system: withdraw the guide catheter and the complete "AVANTGARDE CHRONO™ CARBOSTENT" device together, as if they were a single unit. An unexpanded stent may be retracted into the guide catheter once only. Other movements inside and out of the distal end of the guide catheter should not be continued. Applying excessive force and/or incorrectly handling the system may lead to release of the stent or cause damage to the positioning catheter.
- Implantation of a coronary stent may cause dissection of the vessel distal and/or proximal to the stent and could also cause acute occlusion of the vessel, making a further operation necessary (CABG, further dilation, placement of additional stents or other procedures).
- During catheter inflation, do not exceed the Rated Burst Pressure.
- Never use air or any other gas substances to inflate the balloon.
- If excessive vessel tortuosity associated with proximal atherosclerosis is observed, the operator could find resistance to device feed. In this case, incorrect handling could cause dissection or rupture of the coronary vessel.

PRECAUTIONS

- Be particularly careful not to handle the stent on the balloon in any way whatsoever. This is very important while the stent protection is being removed, while the guidewire is being inserted and while it is being fed through the haemostatic valve.
- The use of an inflation device fitted with a pressure gauge is recommended.
- To avoid damaging the stent, take great care whenever passing the guidewire or the balloon catheter across a stent that has just been expanded.
- When multiple stents are needed, the stent materials must be of similar composition.
- An accurate choice of the stent size is important in order for the procedure to be successful. Generally speaking, the size of the stent must correspond to the diameter of the vessel in question and to the length of the lesion. A slightly larger stent is preferable to one smaller in size.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that CoCr stent product lines are MR Conditional. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5-Tesla (1.5 T) or 3-Tesla (3 T).
- Maximum spatial field gradient of 2.890 G/cm (28.90 T/m)

Under the scan conditions defined below, CoCr stent product lines in a single stent configuration are expected to produce a maximum temperature rise of less than 2.4 °C in a 1.5 T MRI system and 2.9 °C in a 3 T MRI system.

	1.5 T	3 T
MR system reported, whole body averaged SAR	2.9 W/kg	2.9 W/kg
Calorimetry measured values, whole body averaged SAR	2.1 W/kg	2.7 W/kg
Highest temperature change	2.4 °C	2.9 °C

In non-clinical testing, the image artifact caused by the device extends approximately 1.0 cm from CoCr stents when imaged with a gradient-echo pulse sequence in a 3 T MRI system.

DRUG SYSTEM

Studies in the clinical literature indicate the need to administer anticoagulation therapy during procedure and antiaggregant treatment in the period after the procedure itself. Given the complexity of the procedure and the many factors influencing the patient's condition, definition of the appropriate therapy is left to the discretion of the physician.

DIRECTIONS FOR USE

Stent + positioning catheter preparation

Prior to the angioplasty procedure, administer Heparin and check that the patient's Activated Clotting Time (ACT) is greater than 300 s.

The stent should be implanted in the target lesion using the positioning catheter on which it is mounted

- a) After having examined the pack for any damage, remove the "AVANTGARDE CHRO-

NO™ CARBOSTENT™ device in a sterile environment.

- b) Check that the positioning catheter does not have kinks, bends or other damage.
- c) Carefully remove stent cover. Check that the stent is not damaged and that it is still to the centre of the balloon.

Any damage to the “AVANTGARDE CHRONO™ CARBOSTENT™” device may impair its performance. If the stent is out of place or damaged, do not use it.

- d) Flush balloon catheter guidewire lumen with a heparin/saline mixture.
- e) Prepare the inflation device according to the manufacturer's instructions.
- f) Purge the air from the balloon on which the stent is mounted, as follows:
 - 1) Fill the inflation device with 4 ml contrast medium;
 - 2) After connecting the inflation device to the positioning catheter Luer connector, point the distal tip of the latter (balloon) vertically downwards;
 - 3) Apply negative pressure and suction for at least 30 seconds. **Allow the pressure to rise back to neutral**, allowing the contrast liquid to fill the system
 - 4) Without introducing air, repeat step 3 and suction for 10-15 seconds until air bubbles no longer appear.

Stent insertion

5F (inner diameter 1.47 mm) or larger guide catheters are compatible for use with “AVANTGARDE CHRONO™ CARBOSTENT™” devices. Guidewires with diameter of 0.014 inches (0.365 mm) or less are compatible for use with the “AVANTGARDE CHRONO™ CARBOSTENT™” device. The choice of guidewire stiffness and tip configuration will depend on the physician's clinical experience.

- a) Keep the positioning catheter on which the stent mounted at ambient pressure during insertion.
- b) Flush the exposed portion of the guidewire with a heparin/saline mixture to remove traces of blood and contrast liquid from it.
- c) Ensure that the haemostatic valve is completely open before inserting the system.
- d) Slowly feed the system over the guidewire, until the stent reaches the site to be treated.
- e) When feeding the system into the target vessel, make sure, by means of fluoroscopy, that the stent and the balloon move as a single unit, checking the position of the respective radio-opaque markers.

WARNING: If unusual resistance is encountered at any time, do not force the system: withdraw the guide catheter and complete stent - positioning catheter system together, as if they were a single unit. Applying excessive force and/or incorrect handling may lead to release of the stent or cause damage to the positioning catheter.

Stent implant and expansion

- a) Position the stent across the target lesion.
- b) Expand the stent by slowly inflating the balloon to the nominal diameter. Table 2 lists the diameters of the system as the balloon inflation pressure varies (column I) for the seven nominal diameter groups: 2.25 mm (column II), 2.50 mm (column III), 2.75 mm (column IV), 3.00 mm (column V), 3.50 mm (column VI), 4.00 mm (column VII) and 4.50 mm (column VIII). The diameters shown on a dark background refer to pressures above the Rated Burst Pressure*. The stent expands at a balloon inflation pressure of about 6 atm.

Stent elastic recoil is between 2 and 7%* depending on stent model and expansion diameter.

* NOTE: These values are reported from results of in-vitro testing.

- c) Deflate the balloon, then verify angiographically that the stent has fully expanded. If necessary, reinflate the balloon in order to obtain an optimum diameter.
- d) Before withdrawing the positioning catheter, deflate the balloon and maintain negative pressure for at least 30 seconds, following the normal procedures used in PTCA.
- e) If necessary, the stent may subsequently be dilated using a balloon for PTCA. However, extreme care must be taken:

**not to expand the 2.25 mm nominal diameter stent over 2.55 mm
not to expand the 2.50 mm nominal diameter stent over 3.05 mm;
not to expand the 2.75 mm nominal diameter stent over 3.05 mm;
not to expand the 3.00 mm nominal diameter stent over 3.85 mm;
not to expand the 3.50 mm nominal diameter stent over 3.85 mm;
not to expand the 4.00 mm nominal diameter stent over 5.05 mm;
not to expand the 4.50 mm nominal diameter stent over 5.05 mm.**

The final diameter of the stent must be suitable to the reference vessel diameter. **Ensure that the stent is not under expanded.**

SIDE EFFECTS

Stent implant may give rise to following side effects:

- Acute myocardial infarction
- Cardiac arrhythmia including ventricular fibrillation
- Unstable angina pectoris
- Dissection, perforation, rupture of the coronary artery
- Distal and/or stent embolism
- Restenosis
- Artefacts during MRI
- Development of allergies
- Death
- Stent occlusion
- Vessel occlusion
- Arterial spasm
- Infection
- Access site haematoma
- Systemic haemorrhage

LIABILITY AND WARRANTY

The Manufacturer guarantees that this device has been designed, manufactured and

packed with the utmost care, using the techniques considered most appropriate among those available in the current state of technology and applying integrated safety standards in the design and manufacture that will guarantee its safe use under the conditions and for the purposes stipulated, observing the precautions described in the paragraphs above, and that will nevertheless reduce, as far as possible, but not eliminate completely, the risks related to the use of the device. The device must be used only under a specialist physician's responsibility, and taking into account any residual risks or possible side effects and complications of the treatment for which it is intended, among those described in other sections of this instruction booklet. Given the technical complexity, the critical nature of treatment choices and the methods used to apply the device, The Manufacturer cannot be held responsible, either explicitly or implicitly, for the successful result following the use of the device or for its effectiveness in resolving a patient's condition. In fact, the results, in terms of both the patient's clinical status of the patient and the functionality and the working life of the device, depend on many factors beyond the manufacturer's control, among which are the patient's conditions, the surgical implant procedure or use, and the methods of handling the device after the pack has been opened. In the light of these factors, therefore, the Manufacturer is solely responsible for the replacement of any device which, upon delivery, is found to have manufacturing defects. To this end, the purchaser must return the device to the Manufacturer, which reserves the right to inspect the device deemed to be faulty and to determine whether the device actually has manufacturing or material defects. The warranty consists exclusively of the replacement of the device found to be faulty with another device of the same or similar type. The warranty only applies if the device is returned correctly packed to the Manufacturer, and is accompanied by a detailed report in writing describing the defects claimed, and, if the device has been implanted, stating the reasons for its removal from the patient. On replacement of the device, The Manufacturer shall reimburse the purchaser the expenses incurred for replacement of the device acknowledged to be faulty. The Manufacturer does not assume any liability in any case of failure to observe the methods of use and the precautions stated in this instruction booklet and in the event of use of the device after the expiry date stated on the pack. Furthermore, the Manufacturer does not assume any liability in relation to the consequences arising from the medical choices and the methods of use or the application of the device; The Manufacturer shall, therefore, not be liable for any damage of any nature, material, biological or moral, as a result of the application of the device, and the choice of implant technique used by the operator. The agents and the representatives of the Manufacturer are not authorized to amend any of the conditions of this warranty or to undertake any further commitments or to offer any guarantee in relation to this product, beyond the terms stated above.

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BLANK PAGE

BLANK PAGE


TABLE 1		STENT LENGTH (mm)							
		7mm	8mm	12mm	16mm	20mm	24mm	25mm	31mm
NOMINAL STENT DIAMETER (mm)	2.25 mm	ICV9601	-	ICV9602	ICV9603	ICV9604	ICV9605	-	-
	2.5 mm	-	ICV9606	ICV9607	ICV9608	ICV9609	-	ICV9610	-
	2.75 mm	-	ICV9611	ICV9612	ICV9613	ICV9614	-	ICV9615	-
	3.0 mm	-	ICV9616	ICV9617	ICV9618	ICV9619	-	ICV9620	ICV9621
	3.5 mm	-	ICV9622	ICV9623	ICV9624	ICV9625	-	ICV9626	ICV9627
	4.0 mm	-	ICV9628	ICV9629	ICV9630	ICV9631	-	ICV9632	ICV9633
	4.5 mm	-	-	ICV9634	ICV9635	ICV9636	-	ICV9637	ICV9638

TABLE 2		BALLOON DIAMETER (mm)							
Inflation pressure (atm)									
	2.25	2.50	2.75	3.00	3.50	4.00	4.50		
6	2.08	2.31	2.54	2.77	3.26	3.74	4.15		
7	2.13	2.37	2.61	2.85	3.35	3.84	4.27		
8	2.19	2.43	2.68	2.93	3.43	3.92	4.38		
9	2.25	2.50	2.75	3.00	3.50	4.00	4.50		
10	2.31	2.55	2.83	3.08	3.58	4.09	4.61		
11	2.36	2.61	2.90	3.15	3.66	4.17	4.68		
12	2.40	2.65	2.95	3.20	3.71	4.25	4.74		
13	2.44	2.70	3.01	3.25	3.77	4.32	4.80		
14	2.48	2.74	3.06	3.29	3.82	4.37	4.86		
15	2.51	2.77	3.09	3.32	3.86	4.41	4.92		
16	2.55	2.80	3.12	3.36	3.90	4.46	4.98		
17	2.59	2.83	3.15	3.39	3.94	4.50	5.04		
18	2.62	2.86	3.19	3.43	3.98	4.54	5.11		
19	2.65	2.89	3.23	3.47	4.03	4.59	5.19		
20	2.68	2.93	3.27	3.51	4.08	4.65	5.27		
21	2.71	2.96	3.30	3.55	4.12	4.69	5.33		
22	2.74	2.99	3.34	3.60	4.18	4.74	5.40		

Pressure above Rated Burst Pressure
Do not exceed Rated Burst Pressure

These values result from in vitro testing

If postdilatation is required, a PTCA balloon can be used. Use extreme care:
not to expand the 2.25 mm nominal diameter stent over 2.55 mm
not to expand the 2.50 mm and the 2.75 mm nominal diameter stents over 3.05 mm
not to expand the 3.00 mm and the 3.50 mm nominal diameter stents over 3.85 mm
not to expand the 4.00 mm and the 4.50 mm nominal diameter stents over 5.05 mm