Ephesos™ II

BARE METAL STENT

EN	ENGLISH.....	4
TR	TÜRKÇE.....	5
FR	FRANÇAIS.....	6
DE	DEUTSCH.....	8
NL	NEDERLANDS.....	9
ES	ESPAÑOL.....	11
IT	ITALIANO.....	13
PT	PORTUGUÊS.....	14
HR	HRVATSKI.....	16
HU	MAGYAR.....	17
BS	BOSANSKI.....	18
CS	ČEŠTINA.....	20
DA	DANSK.....	21
SR	SRPSKI.....	22
SV	SVENSKA.....	24
RU	РУССКИЙ.....	25
PL	POLSKI.....	27
RO	ROMÂNĂ.....	28
MK	МАКЕДОНСКИ.....	30
BG	БЪЛГАРСКИ.....	31
EL	ΕΛΛΗΝΙΚΑ.....	32
SK	SLOVENČINA.....	34
ZH	简体中文.....	35
JA	日本語.....	36
AR	العربية.....	38
FA	فارسی.....	39

Rx only

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EPHESOS™ II CORONARY STENT DELIVERY SYSTEM

DESCRIPTION

The Alvimedica EPHESOS™ II Stent is a balloon-expandable stent, laser-cut from 316 LVM stainless steel tube. The EPHESOS™ II Stent is available in a range of sizes, namely 6 mm, 8 mm, 9 mm, 12 mm, 15 mm, 18 mm, 20 mm, 25 mm, 28 mm and 32 mm. The EPHESOS™ II Stent is intended for use mounted on a TURQUOISE™ PICA Balloon Dilatation Catheter.

Stent Inner Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter Compatibility ID	Nominal Pressure (bar)	Rated Burst Pressure (bar)
2.00	9-12-15-18-20	5F (1.422 mm)	8	16
2.50	9-12-15-18-20-25-28-32	5F (1.422 mm)	8	16
2.75	9-12-15-18-20-25-28-32	5F (1.422 mm)	8	16
3.00	9-12-15-18-20-25-28-32	5F (1.422 mm)	8	16
3.50	9-12-15-18-20-25-28-32	5F (1.422 mm)	8	16
4.00	6-8-9-12-15-18-20-25-28-32	5F (1.422 mm)	8	16
4.50	6-8-9-12-15-18-20	5F (1.422 mm)	8	16

STENT DELIVERY SYSTEM

The Alvimedica TURQUOISE™ PICA Balloon Dilatation Catheter has been designed to dilate stenotic atherosclerotic lesions in coronary arteries or bypass grafts. The dilation part of the catheter is the balloon near the distal tip. The catheter is hydrophilic-coated from the distal tip towards the proximal end over a length of 105 cm, excluding the effective length of the balloon and 45 cm of the proximal shaft. Radiopaque marker tapes are located on both the proximal and distal shoulders of the balloon. A separate lumen on the catheter shaft is intended for use as a guidewire lumen beginning at approximately 27 cm from the distal tip. The proximal tip of the catheter is intended for use as the balloon inflation port. The balloon is inflated by injecting contrast material from this tip. The balloon material can be inflated to a given size at a given pressure.

The maximum diameter of the guidewire that can be used with the TURQUOISE™ Balloon Dilatation Catheter is 0.014" (0.356 mm).

Warning: EPHESOS™ II Stent reaches its nominal dilated diameter when the delivery system balloon is inflated at a pressure of 8 bars. The delivery system balloon should not be inflated at pressures that dilate the stent beyond the intended vessel lumen diameter or at pressures that exceed the pre-determined maximum pressure.

INDICATIONS

The EPHESOS™ II Stent system is indicated for:

- Patients who have symptomatic ischemic heart disease due to de novo stenosis or restenosis of coronary arteries and who are eligible for PICA procedures
- Treatment of acute obstruction complicating PICA procedures
- Cases of unsuccessful PICA procedures
- Decreasing the restenosis ratio
- Patients that are eligible for coronary bypass surgery

CONTRAINDICATIONS

The EPHESOS™ II Stent system is contraindicated in:

- Patients who cannot receive anti-platelet or anti-coagulant therapy
- Patients who are not eligible for coronary bypass surgery
- Patients who have diffuse distal vessel disease that would decrease blood flow through the stent lumen
- Patients with an unprotected left main coronary artery
- Ostial lesions
- Lesions with coronary embolism
- Patients with a known metal allergy
- Pregnant women
- Tortuous lesions with angles of more than 60 degrees.

WARNING

- This product has been designed for single use only. Do not resterilize

or reuse this device.

- All procedures are to be performed once the catheter has been introduced into the body and should be carried out under high-quality fluoroscopy. Never pull or push the catheter unless the balloon has been fully deflated under vacuum. If any resistance is encountered during the procedure, simply stop and try to identify the cause and then gradually advance the balloon catheter to prevent the occurrence of kinking on the proximal shaft. If you fail to identify the cause, remove the system as a single unit (stent, delivery system and guiding catheter).
- Do not exceed the rated burst pressure as indicated in the Instructions for Use when inflating the balloon. To monitor the pressure, use of an inflation device with an incorporated pressure gauge is recommended.
- Delivery of the stent should only be performed in medical centers capable of carrying out emergency coronary bypass surgery in case of life-threatening complications.
- Do not use the stent and the delivery system after the expiry date indicated on the label.
- Alvimedica is not liable for any damage or injury that may be caused by the reuse of the stent and the delivery system.
- Do not try to remove the stent wedged onto the delivery system and do not make any changes to the procedure. The stent should not be removed from the delivery system and placed in another balloon catheter, as this may damage the balloon or the stent.
- Do not use the contrast materials Ethiodol or Lipiodol.
- Do not expose the delivery system to organic solvents (i.e. alcohol, etc.).
- Do not inflate the balloon at a higher pressure than is required.
- An inappropriate anticoagulant and/or antiplatelet therapy may cause stent thrombosis.
- To prevent the possibility of metal corrosion, do not place stents made of different metals inside one another, next to each other or in direct contact with each other.
- Do not use air or any other gas to inflate the balloon.

The product might be a potential biohazard, after use. The product must be handled and disposed of in accordance with accepted medical practice, applicable local laws and regulations. Alvimedica is not responsible for the handling and disposal of the product after use.

PRECAUTIONS

- Prior to the stent delivery, check whether the shape and size of the stent and the delivery system to be used are suitable for the specific case.
- The stent and the delivery system should only be used by experienced physicians who have been adequately trained in angiography, stent delivery and PICA operations.
- Do not use if the inner packaging has been opened or damaged.
- Do not sterilize in an autoclave.
- Keep in a cool, dry and dark place.
- Prior to the stent delivery, the lesion should be pre-dilated using standard PICA techniques.
- The delivered stent can be damaged by physical contact with auxiliary materials.
- Repeat dilation of the delivered stent may cause dissection. This can be corrected by delivering another stent. In cases where more than one stent is delivered, their tips should be overlapped.
- When a fully or partially dilated stent needs to be introduced, care should be taken to avoid deforming or damaging the stent.
- It is recommended that mechanical atherectomy and the use of laser catheters be avoided in the area where the stent has been delivered.
- The stent may cause artifacts in MR imaging due to the distortion in the magnetic field. Moreover, in order to decrease the risk of stent migration in powerful magnetic fields, MR imaging should not be performed until the stent has been completely covered by endothelial tissue (i.e. a period of about 8 weeks).
- It is recommended that hyperthermia be avoided in the area where the stent has been delivered.

ADVERSE EFFECTS / COMPLICATIONS

The following adverse effects may occur during or subsequent to stent delivery:

- Acute myocardial infarction
- Unstable angina
- Complete obstruction
- Tearing or perforation of, or damage to the coronary artery
- Coronary artery dissection, perforation, injury or tearing
- Re-stenosis of the stent-delivered vessel
- Spasm of the coronary artery
- Deformation of stent symmetry during delivery or while introducing another device (e.g. a catheter, balloon, guidewire, etc.) through the stent
- Hemorrhage or hematoma
- Arrhythmias, including ventricular fibrillation
- Drug reactions, allergic reaction to contrast material
- Hypo-/ hypertension
- Infection
- Arteriovenous fistula
- Embolism
- Urgent coronary artery bypass graft surgery
- Gastrointestinal bleeding due to anticoagulant/antiplatelet therapy
- Intimal tearing
- Formation of pseudo-aneurysm
- Stent migration
- Thrombosis
- Tissue necrosis
- Death

PREPARATION AND TESTING PRIOR TO USE

Ensure that the diameter of the stent to be used is identical to the diameter of the treated vessel.

Auxiliary material (not included in the package)

- Guiding catheter(s) with a 5F diameter (1.422 mm) or larger
- Syringes
- Guidewires 0.014" (0.356 mm) in diameter or less
- Inflation device
- Introducer sheath

Preparation

1. Gently remove the delivery system from the packaging. Be careful when taking the stylet and the protection tube out of the distal tip of the catheter. Check the integrity of the wedged stent and the centering on the balloon; also check whether there are any projections on the connections.

Note: The stent should be centered on the balloon.

Warning: Do not use gauze, as its fibers may cling to and impair the stent.

2. Fill the 20 cc syringe with 10 cc of saline solution. Place a needle on the syringe and carefully insert the needle through the distal tip of the catheter and flush the guidewire lumen.
3. Hold the stent carefully, maintaining the position of the stent wedged on the balloon, while taking the catheter out of the packaging. Carefully load it onto the guidewire, insert it into the guiding catheter and pass it through the Y-connector with the hemostatic valve, paying the utmost attention during each of these stages.
4. Fill the 20 cc syringe with a 50/50% mixture of contrast material and saline solution. Purge the syringe completely of air, connect the syringe to the stopcock and fill the stopcock lumen with the contrast-saline mixture.
5. Prepare the inflation device in accordance with the manufacturer's instructions.
6. Connect the stopcock to the inflation port on the proximal part of the catheter. Connect the inflation device to the other arm of the stopcock. While carrying out these procedures, tighten the connections to prevent air from entering the system.

7. Switch the stopcock for the inflation device to the "off" position.
8. While holding the syringe vertically with the piston positioned above, pull the piston and aspirate until no air bubbles remain in the syringe.
9. Switch the stopcock for the inflation device to the "on" position.

Note: Do not apply negative pressure on the system with the inflation device.

DELIVERY AND DEPLOYMENT PROCEDURES

Note: It is not recommended that the Alvimedica EPHESOS™ II Stent and the delivery system be used for pre-dilation purposes.

1. Dilate the lesion using the PTCA balloon. It is preferable to match the diameters of the balloon and the vessel.
2. Advance the prepared stent and the delivery system over the guidewire (max. diam. 0.014" (0.356 mm)) towards the Y-connector with the hemostatic valve (inner diameter > 0.074" (1.88 mm)). The valve should be completely open while the stent and delivery system are introduced.
- Note:** If you encounter any resistance, do not force passage, as this may damage the stent or cause it to dislodge from the balloon. Ensure that the guiding catheter is in the correct position and that its tip is stationary.
3. Advance the stent inside the guiding catheter under direct fluoroscopic visualization. Using standard PTCA techniques, position the stent across the lesion.
4. Expand the stent in the lesion using standard angioplasty techniques. While holding the balloon in a deflated position, check the expansion of the stent using standard angiographic imaging. An inflation pressure of 8 bars is sufficient to expand the stent.

Warning: Expanding the stent beyond the limit and exerting more inflation pressure than is recommended may lead to dissection.

REMOVAL OF THE STENT DELIVERY SYSTEM

1. After the stent has been delivered, exert negative pressure and deflate the balloon. Once the balloon has been completely deflated, carefully withdraw the delivery system from the guiding catheter under fluoroscopy while maintaining the guidewire position.
2. Repeat angiographic imaging to assess the stented area. If the stent is under-dilated, carefully position the original delivery system or another catheter with appropriate balloon diameter in the stent and repeat the previous steps to expand the stent to the required diameter.

TERMS OF THE WARRANTY

Alvimedica guarantees that each and every component of this product has been manufactured, packed, tested and sterilized without any defect to its workmanship or material. Each product has been tested prior to packaging. Alvimedica shall exchange any product with manufacturing or packaging defect(s) with a new one upon the decision of a competent binding judicial authority. Due to biological variations among individuals, no product is 100 % effective in every case. Therefore, and since Alvimedica has no control over conditions under which the product is used after sales, the selection of patients and the methods of application; Alvimedica disclaims all warranties with respect to the product, expressed or implied including, without limitation, those of merchantability, and fitness for a particular purpose. Alvimedica is not directly or indirectly responsible for any injury or damage to or loss of any person resulting from the use of the product, nor is Alvimedica responsible or liable, directly or indirectly, for any injury, damage or loss that may result from re-use or re-sterilization.

LIMITATION OF LIABILITY

In no event shall Alvimedica be liable for any special, consequential, indirect, collateral, incidental or punitive damages or lost profits or failure to realize expected savings or other commercial or economic loss of any kind, arising out of any breach of this IFU including breach of any precautions, warnings and/or any recommendation mentioned herein.

INTELLECTUAL PROPERTY PROTECTION

All trademarks, logos, inventions, know-how, technology, proprietary information and other intellectual property rights of this product belong solely to Alvimedica. Alvimedica reserves all the rights to demand, sue, seek relief before judicial authorities in case of an infringement regarding any unauthorized use or misuse of intellectual property or reverse engineering related to the product.

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