


# Commander™

CoCr BARE METAL STENT

Very thin cobalt chromium coronary stent  
designed to provide homogenous vessel coverage


- Flexible Cobalt Chromium alloy
- Homogenous vessel coverage
- 70µm thin stent strut




# Commander™

CoCr BARE METAL STENT


## UNIQUE FEATURES



Flexible Cobalt Chromium alloy




Homogenous vessel coverage



70µm thin stent strut

## TECHNICAL DATA



CATHETER SPECIFICATIONS	
Category:	Rapid Exchange balloon delivery system
Guiding catheter compatibility:	5F
Guidewire compatibility:	0.014"
Usable catheter length:	150 cm
Proximal shaft:	stainless steel AISI 304, hypotube with PEBA jacket
Distal shaft:	trilayered polyamide blend shaft with reinforced inner lumen and Hydrophilic coating

STENT SPECIFICATIONS	
Category:	BMS
Material:	L605 (CoCr alloy)
Strut thickness	70 µm
X-ray visibility:	Moderate and homogeneous
Foreshortening	< 2%

BALLOON SPECIFICATIONS	
Nominal Pressure:	8 atm
Rated Burst Pressure:	16 atm

## ORDER INFORMATION COMMANDER™

L [mm] \ Ø [mm]	9	12	15	18	20	22	25	28	32
2.00	508110020801	508110020802	508110020803	508110020804	508110020805	508110020806			
2.50	508110020807	508110020808	508110020809	508110020810	508110020811	508110020812	508110020813	508110020814	508110020815
2.75	508110020816	508110020817	508110020818	508110020819	508110020820	508110020821	508110020822	508110020823	508110020824
3.00	508110020825	508110020826	508110020827	508110020828	508110020829	508110020830	508110020831	508110020832	508110020833
3.50	508110020834	508110020835	508110020836	508110020837	508110020838	508110020839	508110020840	508110020841	508110020842
4.00	508110020843	508110020844	508110020845	508110020846	508110020847	508110020848	508110020849	508110020850	508110020851
4.50	508110020852	508110020853	508110020854	508110020855	508110020856	508110020857			

[www.alvimedica.com](http://www.alvimedica.com)

Please contact your local sales representative for product availability. The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2020, Alvimedica. All Rights Reserved.

