


# Ephesos™ II

## BARE METAL STENT

Very thin stainless steel coronary stent designed to provide homogenous vessel coverage

- Flexible 316L SS platform
- Homogenous vessel coverage
- 80µm thin stent strut




# Ephesos™ II

## BARE METAL STENT

### UNIQUE FEATURES



Flexible 316L SS platform




Homogenous vessel coverage



80µm thin stent strut

### TECHNICAL DATA



#### CATHETER SPECIFICATIONS

Category:	Rapid Exchange balloon delivery system
Guiding catheter compatibility:	5F
Guidewire compatibility:	0.014"
Usable catheter length:	150 cm
Proximal shaft:	stainless steel AISI 304, hypotube with PEBA jacket
Distal shaft:	trilayered polyamide blend shaft with reinforced inner lumen and Hydrophilic coating

#### STENT SPECIFICATIONS

Category:	BMS
Material:	316L stainless steel
Strut thickness:	80 µm
X-ray visibility:	Moderate and homogeneous
Foreshortening:	< 2%

#### BALLOON SPECIFICATIONS

Nominal Pressure:	8 atm
Rated Burst Pressure:	16 atm

### ORDER INFORMATION EPHEOS™ II

∅ [mm]	6	8	9	12	15	18	20	25	28	32
2.00			508110060301	508110060302	508110060303	508110060304	508110060305			
2.50			508110060306	508110700607	508110060308	508110060309	508110060310	508110060311	508110060312	508110060313
2.75			508110060347	508110060348	508110060349	508110060350	508110060351	508110060352	508110060353	508110060354
3.00			508110060314	508110060315	508110060316	508110060317	508110060318	508110060319	508110060320	508110060321
3.50			508110060322	508110060323	508110060324	508110060325	508110060326	508110060327	508110060328	508110060329
4.00	508110060330	508110060331	508110060332	508110060333	508110060334	508110060335	508110060336	508110060337	508110060338	508110060339
4.50	508110060340	508110060341	508110060342	508110060343	508110060344	508110060345	508110060346			

[www.alvimedica.com](http://www.alvimedica.com)

Please contact your local sales representative for product availability. The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2020, Alvimedica. All Rights Reserved.

