


Fluydo™ NC

NON COMPLIANT PTCA RX BALLOON DILATATION CATHETER

Fluydo™ NC is a non compliant balloon designed to deliver the strength you need where it's needed


- Versatile non-compliant balloon with superb crossability
- Perfect combination of high RBP with ultra low compliance
- 6F kissing balloon technique fully compatible*



Fluydo™ NC

NON COMPLIANT PTCA RX BALLOON DILATATION CATHETER

TECHNICAL DATA



CATHETER SPECIFICATIONS		BALLOON SPECIFICATIONS	
Category:	Rapid Exchange balloon catheter	Category:	Non Compliant
Guiding catheter compatibility:	5F	Material:	Polyamide
Guidewire compatibility:	0.014"	X-ray visibility:	Two platinum markers
Usable catheter length:	142 cm	Nominal Pressure:	12 atm
Proximal shaft:	stainless steel AISI 304, low friction hypotube with Teflon coating	Rated Burst Pressure:	22 atm (1.75mm to 4.00mm) 20 atm (4.50mm) 19 atm (5.00mm)
Distal shaft:	mix Polyamide/PEBA with Hydrophilic coating	Average Burst Pressure:	28 atm (1.75mm to 4.00mm) 26 atm (4.50mm to 5.00mm)
Tip material /entry profile:	PEBA blend / 0.017"		

COMPLIANCE CHART

Pressure [atm] \ Ø	1.75 mm	2.00 mm	2.25 mm	2.50 mm	2.75 mm	3.00 mm	3.25 mm	3.50 mm	3.75 mm	4.00 mm	4.50 mm	5.00 mm
12	1.75	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
14	1.77	2.03	2.28	2.53	2.79	3.04	3.29	3.55	3.80	4.05	4.56	5.07
16	1.80	2.05	2.31	2.56	2.82	3.08	3.33	3.59	3.85	4.10	4.62	5.13
18	1.82	2.08	2.34	2.60	2.86	3.11	3.38	3.64	3.89	4.16	4.68	5.20
19	1.83	2.09	2.35	2.61	2.88	3.13	3.40	3.66	3.92	4.18	4.71	5.23
20	1.85	2.10	2.37	2.63	2.89	3.15	3.42	3.68	3.94	4.21	4.74	
22	1.87	2.13	2.40	2.66	2.93	3.19	3.46	3.73	3.99	4.26		

Nominal Pressure

Rated Burst Pressure

ORDER INFORMATION FLUYDO™ NC

L [mm] \ Ø [mm]	6	8	12	15	20	30
1.75	ICNC17506	ICNC17508	ICNC17512	ICNC17515	ICNC17520	-
2.00	ICNC20006	ICNC20008	ICNC20012	ICNC20015	ICNC20020	-
2.25	ICNC22506	ICNC22508	ICNC22512	ICNC22515	ICNC22520	-
2.50	ICNC25006	ICNC25008	ICNC25012	ICNC25015	ICNC25020	ICNC25030
2.75	ICNC27506	ICNC27508	ICNC27512	ICNC27515	ICNC27520	-
3.00	ICNC30006	ICNC30008	ICNC30012	ICNC30015	ICNC30020	ICNC30030
3.25	ICNC32506	ICNC32508	ICNC32512	ICNC32515	ICNC32520	-
3.50	ICNC35006	ICNC35008	ICNC35012	ICNC35015	ICNC35020	ICNC35030
3.75	ICNC37506	ICNC37508	ICNC37512	ICNC37515	ICNC37520	-
4.00	ICNC40006	ICNC40008	ICNC40012	ICNC40015	ICNC40020	-
4.50	-	ICNC45008	ICNC45012	ICNC45015	ICNC45020	-
5.00	-	ICNC50008	ICNC50012	ICNC50015	ICNC50020	-

Manufactured by: CID S.p.A.
member of Alvimedica group
Via Crescentino, sn - 13040 Saluggia (VC), Italy

*6F kissing balloon technique fully compatible: any two balloons with at least one with diameter ≤ 3,25mm

www.alvimedica.com

Please contact your local sales representative for product availability. The images shown within this document are of a purely illustrative nature and are not intended to depict reality.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions and instruction for use.

©2020, Alvimedica. All Rights Reserved.

